# Tesmed+

## TE-880R PLUS

MANUALE D'USO - USER MANUAL  
NOTICE D'INSTRUCTIONS  
BEDIENUNGSANLEITUNG - MANUAL DE USO

Per ottenere una copia del manuale in formato elettronico  
To obtain an electronic copy of the manual  
Pour obtenir une copie électronique du manuel  
So erhalten Sie eine elektronische Kopie des Handbuchs  
Para obtener una copia electrónica del manual

<http://www.tesmed.it/pages/download-manual>  
or write to  
[info@tesmed.com](mailto:info@tesmed.com)

ITALIANO

ENGLISH

FRANÇAIS

DEUTSCH

ESPAÑOL

ITALIANO

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Tesmed does not guarantee its contents and reserves the right to improve and amend it at any time without prior notice. However, Amendments will be published in a new edition of this manual.

All Rights Reserved.

Tesmed TE-880R PLUS Rev.V1.1 © 2024,  
revised on August 2024.

**Declaration of conformity:**

Shenzhen Roundwhale Technology Co., Ltd. declares that the device complies with following normative documents:

IEC60601-1, IEC60601-1-2, IEC60601-1-11, IEC60601-2-10,  
IEC62304, ISO10993-5, ISO10993-10, ISO10993-23,  
ISO10993-1, ISO14971

## TABLE OF CONTENTS

1. FOREWORD.....	49
2. SAFETY INFORMATION.....	51
(important safety information, contraindications and warnings)	
3. GETTING TO KNOW YOUR DEVICE.....	56
(accessories, display and device illustration, description of the function of each individual button)	
4. TECHNICAL INFORMATION.....	59
5. OPERATING INSTRUCTION.....	60
(electrode and cable connection, electrode positioning)	
6. INSTRUCTIONS FOR USE .....	66
7. RECOMMENDATIONS FOR USE .....	74
8. PROGRAMS AND RELATED INDICATIONS .....	75
9. CLEANING AND MAINTENANCE.....	79
10. TROUBLESHOOTING.....	80
11. STORAGE.....	82
12. DISPOSAL.....	83
13. ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES .....	84
14. NORMALIZED SYMBOLS.....	90
15. WARRANTY .....	91

TESMED copyright - It is absolutely forbidden the partial or total reproduction of any present content

Be sure to read this instruction manual before operating and keep it where safe.

## 1. FOREWORD

### Introduction

The device Tesmed TE-880R PLUS is a dual channel output TENS, EMS and MASSAGE stimulator. Before using, please read all the instructions in this user manual carefully and keep it safe for future use.

The TE-880R PLUS stimulator belongs to the group of electrical stimulation systems. It has three basic functions- TENS (Transcutaneous Electrical Nerve Stimulation), EMS (Electronic Muscle Stimulation) and MASSAGE.

Function of the TE-880R PLUS stimulator: The device has 73 programs (11 TENS programs, 10 EMS programs and 52 MASSAGE programs) and applies electric currents in the lowfrequency range for therapy. Each program controls the generated electric impulses, their intensity, frequency and pulse width.

Based on simulating the body's natural pulses, the mechanism of electrical stimulation equipment is to create electric impulses that are transcutaneous transmitted to nerves or muscle fibers through the electrode. The intensity of the dual channel can be adjusted independently and applied individually to one body part. This dual channel device can be used with eight pieces of electrodes. The electrical pulse is firstly transmitted to the issue, then it affects the transmission of stimulation in nerves as well as muscle tissues in the body parts.

## **1.2 Medical background**

### **1.2.1 ABOUT PAIN**

Pain is an important signal in the human body warning system. It reminds us that something is wrong, without which, abnormal conditions may go undetected, causing damage or injury to vital parts of our bodies.

Aside from its function in diagnosis, longlasting persistent pain serves useless purpose.

Pain does not occur until encoded message travels to the brain where it is decoded, analyzed, and reacted to, from the injured area along the small nerves leading to the spinal cord. There the message is transmitted to different nerves that travel up the spinal cord to the brain. Then the pain message is interpreted, referred to and pain is felt.

### **1.2.2 WHAT IS TENS ?**

TENS (Transcutaneous Electrical Nerve Stimulation) is effective in relief of pain. It is daily used and clinically proven by physiotherapists, caregivers and top athletes around the world. High-frequency TENS currents activates the pain inhibiting mechanisms of the nervous system. Electrical impulses from electrodes, placed on the skin over or near the pain area, stimulate the nerves to block the pain signals to the brain, causing the pain go unperceived. Low-frequency TENS currents facilitate the release of endorphins, the body's natural painkillers.

### **1.2.3 WHAT IS EMS ?**

Electrical Muscle Stimulation is an internationally accepted and proven way of treating muscular injuries. It works by sending electronic pulses to the muscle needing treatment that causes the muscle to exercise passively. It is a product deriving from the square waveform, originally invented by John Faraday in 1831. Through the square wave pattern it is able to work directly on muscle motor neurons. The EMS System has low frequency and this in conjunction with the square wave pattern allows direct work on muscle groupings.

### **1.2.4 WHAT IS MASSAGE ?**

The massage function is non-medical function. The Massage stimulation program provides relaxing muscle vibration to loosen tight muscles.

Trained copyright - It is absolutely forbidden the partial or total reproduction of any present content.

## 2. SAFETY INFORMATION

### 2.1 Intended use

#### Intended purpose

The device is designed to be used for temporary relief of pain, including the acute and chronic pain relief.

The device is designed to be used for stimulate muscles for prevention of muscular atrophy, muscle strengthening and improve, Increasing local blood circulation and facilitate muscle performance.

#### Target population

The device using the object (patient) must be 18 years or older of adults.

#### Intended user

Medical staff or lay persons.

#### Intended condition

Intended for use in the home, hospital and health care facilities.

#### Indications

- 1) It is used for temporary relief of pain associated with sore and aching muscles in the neck, shoulder, back, joint, hip, hand, abdomen, upper extremities (arm) and lower extremities (leg) due to strain from exercise or normal household work activities.
- 2) Relaxation of muscle spasms.
- 3) Prevention or retardation of disuse atrophy.
- 4) Increasing local blood circulation.
- 5) Muscle re-education.
- 6) Immediate post-surgical stimulation of calf muscles to prevent venous thrombosis.
- 7) Maintaining or increasing range of motion.

### 2.2 Important Safety Precautions and Warnings

It is important that you read all the warnings and precautions included in this manual because they are intended to keep you safe, prevent risk of injury and avoid a situation that could result in damage to the device.

## SAFETY SYMBOLS USED IN THIS MANUAL

### 2.2.1 Contraindication

- 1) Do not use this device if you have a cardiac pacemaker, implanted defibrillator, or other implanted metallic or electronic device. Such use could cause electric shock, burns, electrical interference, or death. 
- 2) The device should not be used when cancerous lesions or other lesions are present in the treatment area.
- 3) Stimulation should not be applied over open wounds or rashes, or over swollen, red, infected, or inflamed areas or skin eruptions (e.g. phlebitis, thrombophlebitis, varicose veins, arteriosclerosis obliterans etc.).
- 4) Electrode placements must be avoided that apply current to the carotid sinus region (anterior neck) or transcranially (through the head).  
- 5) Apprehensive patients-usage of Electrical stimulation requires patient cooperation, hence the procedure shouldn't be attempted in patients with a communication handicap or a mental disability.
- 6) Patients with cerebrovascular problems-patients with a history of aneurysm, stroke and transient ischaemia shouldn't be treated using electrical stimulation, as it stimulates peripheral blood flow which can be fatal in such cases.
- 7) Epileptic patients-Electrical stimulation "pulses" have the potential to trigger a seizure.
- 8) Acute pain cases/pain of unknown etiology-usage of TENS in undiagnosed cases may hinder in the diagnosis.
- 9) Do not use in pregnancy, especially in the first trimester.

**2.2.2**  **WARNING**

- 1) If you have had medical or physical treatment for your pain, consult with your physician before use.
- 2) If your pain is not subdued, which becomes more than mild, or lasts for more than five days, stop using the device and consult with your physician.
- 3) Do not apply stimulation over your neck because this could cause severe muscle spasms resulting in closure of your airway, difficulty in breathing, or adverse effects on heart rhythm or blood pressure.
- 4) Do not apply stimulation across your chest because the introduction of electrical current into the chest may cause rhythm disturbances to your heart, which could be lethal.
- 5) Do not apply stimulation in the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms), which may not operate properly when electrical stimulation device is in use.
- 6) Do not apply stimulation when in bath or shower.
- 7) Do not apply stimulation while sleeping.
- 8) Do not apply stimulation while driving, operating machinery, or during any activity when electrical stimulation can put you at risk of injury.
- 9) Apply stimulation only to normal, intact, clean, healthy skin.
- 10) The long-term effects of electrical stimulation are unknown. Electrical stimulation device cannot replace drugs.
- 11) Stimulation should not take place while the user is connected to high-frequency surgical equipment, which may cause burn injuries on the skin under the electrodes, as well as problems with the stimulator.

- 13) Do not use the stimulator in the vicinity of shortwave or microwave therapy equipment, since this may affect the output power of the stimulator.
- 14) Never use it near the cardiac area. Stimulation electrodes should never be placed anywhere on the front of the thorax (marked by ribs and breastbone), but above all not on the two large pectoral muscles. There it can increase the risk of ventricular fibrillation and lead to cardiac arrest. 
- 15) Never use it on the eye, head and face area. 
- 16) Never use it near the genitals.
- 17) Never use it on the areas of the skin which lack normal sensation.
- 18) Keep electrodes separated during treatment. It could result in improper stimulation or skin burns if electrodes are in contact with each other.
- 19) Keep the stimulator out of reach of children.
- 20) Consult your doctor if you are in any doubt whatsoever.
- 21) Discontinue it and do not increase the intensity level if you feel discomfort during use.

### 2.2.3 Precautions

- 1) TENS is not effective for pain of central origin including headache.
- 2) TENS is not a substitute for pain medications and other pain management therapies.
- 3) TENS is a symptomatic treatment and, as such, suppresses the sensation of pain that would otherwise serve as a protective mechanism.
- 4) Effectiveness is highly dependent upon patient selection by a practitioner qualified in the management of pain patients.

- 5) Since the effects of stimulation of the brain are unknown, stimulation should not be applied across your head, and electrodes should not be placed on opposite sides of your head.
- 6) The safety of electrical stimulation during pregnancy has not been established.
- 7) You may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium (silica gel).
- 8) If you have suspected or diagnosed heart disease or epilepsy, you should follow precautions recommended by your physician.
- 9) Caution if you have a tendency to bleed internally, e.g. following an injury of fracture.
- 10) Consult with your physician prior to use the device after a recent surgical procedure, because stimulation may disrupt the healing process.
- 11) Caution if stimulation is intended to be applied over the menstruation or pregnant uterus.
- 12) For single patient use only.
- 13) This stimulator should not be used by patients who is non-compliant and emotionally disturbed including whom with dementia or low IQ.
- 14) The instruction of use is listed and should be obeyed; any improper use may be dangerous.
- 15) Rare cases of skin irritation may occur at the site of the electrode placement following long-term application.
- 16) Do not use this device in the presence of other equipment which sends electrical pulses to your body.
- 17) Do not use sharp objects such as a pencil or ballpoint tip to operate the buttons on the control panel.
- 18) Check the electrode connections before each use.
- 19) Electrical stimulators should be used only with the electrodes recommended for use by the manufacturer.

- 20) When output of device more than 10mA or 10 V, the intensity of Channel will flicker.
- 21) Users should consult a healthcare professional before using the device.
- 22) The user shall report any serious incident related to the device to the manufacturer and the competent authorities of the Member States establishing the user and / or the patient.

### 2.2.4 Adverse Reactions


- 1) Possible skin irritation or electrode burn under the electrodes may occur.
- 2) Possible allergic skin reaction to tape or gel may occur.
- 3) On very rare occasions, first-time users of EMS report feeling light-headed or faint. We recommend that you use the product while seated until you become accustomed to the sensation.
- 4) If symptoms of tachycardia and extrasystolia (rapid heartbeat or extra stimulation) appear during treatment, stop the treatment and seek medical attention immediately.
- 5) If the stimulation makes you uncomfortable, reduce the stimulation Intensity to a comfortable level and contact your physician if problems continue.

## 3. GETTING TO KNOW YOUR DEVICE

### 3.1 Package includes

Description	QTY
The TE-880R PLUS Stimulator	1
Electrode pad (50mmX50mm)	8
Four-core electrode cables	2
Covers for cable plugs	4
User manual	1
Belt clip	1

### 3.2 LCD DISPLAY

No.	Function description	No.	Function description
1	Body part	5	Key locking symbol
2	intensity for Channel 2	6	Program NO.
3	intensity for Channel 1	7	Low battery indicator
4	Treatment time in standby mode or during treatment. Pulse width or pulse rate in the Setting mode	8	Massage symbols
		9	Treatment mode

### 3.3 DEVICE ILLUSTRATION (view pagina 57)

© 2008 TESMED - All rights reserved. TESMED is a registered trademark of TESMED.

No.	Description
1	<b>LCD display</b>
2	<p><b>[ON/OFF/M] button:</b>            At power saving mode (unlit display) press the [ON/OFF/M] button to turn on device.            At standby mode, press the [ON/OFF/M] button to select the treatment mode (TENS, EMS o MASS); at standby mode, press and hold the [ON/OFF/M] button to turn off the device.            At treating mode,(when the intensity has already been selected) press the [ON/OFF/M] button to stop the treatment.            At setting mode, (after selecting the programme) press the [ON/OFF/M] button to enter the standby mode.</p>
3	<p><b>[T] button:</b>            At standby mode, press the [T] button to increase the treatment time (every 5 minutes); Press and hold the [T] button to set the frequency, pulse width and time (for U programmes only)            At setting mode, press and hold [ T] button to select pulse rate, pulse width or treatment time - note : adjust the parameters each time by pressing the + or - key, confirm with the T key, then Press and hold the [T] button to set the frequency, pulse width and time (for U programmes only)</p>
4	<p><b>[CH] button:</b>            At standby mode and treating mode, press the [CH] button to select the treatment channel.</p>
5	<p><b>[-] button:</b>            At treating mode, press the [-] button to decrease the intensity of CH1 and CH2, CH1 or CH2.            At setting mode, press the [-] button to decrease the corresponding data for the pulse rate, pulse width and treatment time.</p>
6	<p><b>[+] button:</b>            At standby or treating mode, press the [+] button to increase the intensity of CH1 and CH2, CH1 or CH2 ;            At setting mode, press the [+] button to increase the corresponding data for the pulse rate, pulse width and treatment time.</p>
7	<p><b>[P] button:</b>            At standby mode, press the [P] button to select the treatment program.            In standby mode, for MASS programmes, press and hold the [P] button to select the desired massage type, each time the button is held down the massage type changes.            At treating mode, press and hold the [P] button turn on/off lock function.</p>
8	<b>Belt Clip</b>
9	<p><b>USB socket and Charger indicator:</b>            When the device is charging, the indicator light will be red            When charging is completed, the indicator light will be green</p>
10	<b>Output socket to which the two cables are connected</b>

TESMED copyright - It is absolutely forbidden the partial or total reproduction of any present content


#### 4. TECHNICAL INFORMATION

Device name	Tesmed TE-880R PLUS
Power sources	3.7 V Li-ion battery
Power supply	Input: 100-240V AC, 50/60Hz, 0.2A; Output: 5V DC, 300mA
Output channel	2
Waveform	Bi-phase square-wave pulse
Output current	Max. 120mA (at 500ohm load)
Output intensity	0 to 40 levels, adjustable
Treatment mode:	TENS, EMS and MASSAGE mode
Number of programs	73 programs TENS: 11 programs EMS: 10 programs MASSAGE: 52 programs
Pulse rate	2Hz ~ 120Hz
Pulse width	50µs ~ 300µs
Treatment time	5 minutes ~ 90 minutes
Operating condition	5°C to 40°C with a relative humidity of 15%-93%, atmospheric pressure from 700 hPa to 1060 hPa
Storage condition	10°C to 55°C with a relative humidity of 10%-95%, atmospheric pressure from 700 hPa to 1060 hPa
Dimension	109*54.5*23mm (L x W x T)
Weight	About 82g
Automatic shutoff	1 minutes
Classification	BF type applied part, internal power equipment, IP22
Size of electrodes pad	50x50mm, square
Output precision	±20% error is allowed for all the output parameters
Device service life	5 years
Electrode pads shelf life	3 years
Software version	V1.0

## 5. OPERATING INSTRUCTION

### 5.1 Connect electrode pads to electrode wires

Insert the electrode wires connector into electrode connector. Make sure they are properly connected to ensure the good performance. Please refer to the picture.


#### **Caution**

Always use the electrode pads which comply with the requirements of the IEC/EN60601-1, ISO10993-1/-5/-10 and IEC/EN60601-1-2, as well as CE and FDA 510(K) regulation.

### 5.2 Connect electrode wires to device

Before proceeding to this step, ensure that the device is completely switched OFF.

Hold the insulated portion of the electrode wire connector, and insert the plug into the receptacle on the top of the main device. Ensure the electrode wires are inserted correctly. The device has two output receptacles controlled by Channel A and Channel B at the top of the unit. You may choose to use one channel with one pair of electrode wires or both channels with two pairs of electrode wires. Using both channels gives the user the advantage of stimulating two different areas at the same time.

### **Caution**

Do not insert the plug of the electrode wires into any AC power supply socket.

## **5.3 Electrode**

### **5.3.1 Electrode options**

The electrodes should be routinely replaced when they start to lose their adhesiveness. If you are unsure of your electrode adhesive properties, please order new replacement electrodes. Replacing electrodes should be reordered under the advice of your physician or the device manufacturer to ensure proper quality. Follow application procedures outlined on electrode packing when using the new replacement electrodes to maintain optimal stimulation and to prevent skin irritation.

### **5.3.2 Place electrodes on skin**

Place the electrode on the body part in need of treatment, according to the instruction of this user manual. Please make the skin clean before use and ensure the skin and electrode connect well.

 **Caution**

1. Always remove the electrodes from the skin with a moderate pull in order to avoid injury in the event of highly sensitive skin.
2. Before applying the self-adhesive electrodes, it is recommended to wash and degrease the skin, and then dry it.
3. Do not turn on the device when the self-adhesive electrodes are not positioned on the body.
4. To remove or move the electrodes, switch off the device or the appropriate channel first in order to avoid unwanted irritation.
5. It is recommended that, at minimum, 40x40 mm self-adhesive square electrodes are used at the treatment area.
6. Never remove the self-adhesive electrodes from the skin while the device is still on.

### 5.4.3 Electrode placement

Tesmed TE-830R PLUS is a kind of OTC stimulator, suitable for home use. You only have to use according to the user manual, place the electrode on the position where you feel pain. Conduct exercise, treatment and adjustment based on your own feeling.

#### Electrode positioning for TENS programmes

##### Attention


The electrodes are white in colour with a green logo. They are shown here with black and red squares that indicate the colour of the pin to which the electrode must be connected.

The squares   shown in the illustrations on the following pages always include two electrodes to which a pair of red and black pins must be connected. Both cables end with two pairs of pins; each pair is comprised of a black and a red pin.

It is also possible to use just two electrodes per cable (instead of four); in this case, connect the two electrodes to one pair of red and black pins, covering the other pair of red and black pins on the cable with the dedicated caps provided.

tesmed copyright - it is absolutely verboten to publish a total reproduction of any part of it

Foot			
Joint (knee)			
Joint (elbow)			
Joint (ankle)			
Joint (wrist)			

#### Position of electrode placement under EMS programs

It is also possible to use just two electrodes per cable (instead of four); in this case, connect the two electrodes to one pair of red and black pins, covering the other pair of red and black pins on the cable with the dedicated caps provided.

Neck		
Shoulder		

Revised copyright - it is absolutely forbidden the partial or total reproduction of any present content


Arm	
Hand	
Back	
Abdomen	
Hip	
Leg	
Foot	

Image copyright. This document prohibits the part or a total reproduction of any content.

## 6. INSTRUCTIONS FOR USE

### 6.1 Turn on

Press the [ON/OFF/M] button to switch on the device; the LCD display lights up. The backlight without operation remains for about 10 seconds after which the display switches to power-saving mode. Pressing the [ON / OFF / M] button switches to standby mode as shown in the following picture:

Attention, the device switches off after about 1 minute of inactivity. **After the first use, if the device is switched on again, it may show the previously selected programme.**


### 6.2 Select treatment mode

Press the [ON/OFF/M] button to select which treatment mode (TENS - EMS - MASS) you will use. The LCD displays as follows:

### 6.3 Select treatment program

Based on your need, press [P] button to select the treatment program. The LCD displays as follows:


Starting from program 2, the area that can be treated will be shown on the display for each type of treatment: a dot will flash on the human figure on the display to indicate the area that can be treated with the selected program.

For example, for the 2 TENS program, a dot will flash on the human figure near the shoulder, to indicate the treatable area. For the list of programs and relative indications see page 75 of the manual.

#### 6.4 For TENS and EMS programmes the default treatment duration is 30 minutes, however, it is possible to adjust the treatment duration for TENS and EMS programmes

The programmes have a default duration of 30 minutes, if you wish to customise the treatment duration, after selecting the desired programme, before giving intensity, press the [T] button (do not hold) to re-adjust the treatment duration. The LCD display shows the following:

#### 6.5 Select treatment channel

By default the two channels are selected, two small arrows appear on the display one to the left and one to the right, with the

[+] button you increase the intensity in the same way on both channels.


Press the [CH] button to select the treatment channel. Channel selection is useful in case you want to set a different intensity level on the two channels. The LCD displays as follows:

#### 6.6 Starting treatment by increasing the intensity

Ensure that the electrodes are placed on the skin and correctly connected.

Press the [+] button to increase the intensity of the selected treatment channel. The LCD displays as follows:


By default both channels are selected, pressing the [+] button will increase the intensity on both channels.


#### 6.7 Adjust the output intensity

Remember to give intensity on both channels.

Press the [+] button to increase output intensity. It will be increased to a higher level after each press. The device has totally 40 levels of output intensity. Please adjust the intensity to the condition that you feel comfortable. The level of output intensity will be shown on the LCD.

At treating status, press and hold [P] button to turn on lock function. The indicator "LOCK" will display on the LCD. This is a safety feature to prevent accidental changes to your settings and to prevent accidentally increasing the output intensity level. Press and hold [P] button to unlock.


**Caution:**

If you feel or become uncomfortable, reduce the stimulation intensity to a more comfortable level and consult with your medical practitioner if problems insist.


**6.8 Stop the treatment and turn off the device**

Press the [ON / OFF / M] button to stop the treatment during treatment mode. Press and hold the [ON / OFF / M] button for 3 seconds if you want to switch off the electro-stimulator and the LCD will switch off as shown below:


#### 6.9 For MASS programmes - how to select the type of massage

In MASS mode, you can select between 4 different massage types according to your preferences. To select between the 4 massage types, press the [P] button for 3 seconds each time. The LCD display shows the following:

#### 6.10 Customisable U programmes - how to set parameters

For the time parameter, simply press the T button to increase the time by 5 minutes at a time. Press and hold the T button to enter the setting mode for pulse frequency, pulse width and treatment time. Once you have selected the parameter you wish to change, press the [+] or [-] button to select the desired value.

Once the desired parameters have been determined, press the ON button and give intensity to start treatment.  
Press [ON/OFF/M] button to return to the standby mode.

#### 6.11 Load detection

If the electrodes do not adhere perfectly to the skin, perhaps because they are worn, or if the leads are not connected to the device, above intensity level 5, the device will detect the anomaly by resetting the intensity and at the same time the symbol, ① or ②, of the reference channel will flash for about 4 seconds. Then the stimulator will return to standby mode as shown below:

#### 6.12 Low battery detection

When the battery is low, the  icon will twinkle to indicate it, stop the device and charge the battery.

#### CHARGING THE BATTERY:

Proceed as follows to recharge the battery:

- This device cannot be used while charging.
- Make sure that the device is no longer connected to the patient (the output cables and electrodes must be disconnected).
- Connect the micro USB cable to the charging port on the device.
- Connect the USB cable to the charger.
- When the device is charging, the indicator light will be red.
- It could take up to 2 hours to reach a full charge.
- When charging is completed, the indicator light will be green.

The life of a rechargeable battery depends on the number of recharging/turnover cycles it undergoes and how these cycles are performed.

The following suggestions will help prolong the life of the battery:

- Whenever the device is not used frequently, charge the battery.
- For longer battery life, discharge the battery as much as possible.

### 6.13 Usage of electrode pads

1. The electrode may only be connected with the COMBO stimulator. Make sure that the device is turned off when attaching or removing the electrode pads.
2. If you want to reposition the electrode during the application, turn the device off first.
3. The usage of electrode may lead to skin irritations. If you experience such skin irritations, e.g. redness, blistering or itching, discontinue using them. Do not use the COMBO stimulator permanently on the same body part, as this may also lead to skin irritations.
4. Electrode pads are private and intended for single person use. Please avoid using them by different persons.
5. The electrode must connect entirely to the skin surface to prevent hot spots, which may lead to skin burns.
6. It is not advisable to use electrodes when they no longer adhere well to the skin. With proper use, an average of about 30 treatments can be made.
7. The adhesive force of the electrodes depends on the skin properties, storage condition, and the number of applications. If your electrode pads no longer fully stick to the skin's surface, replace them with new ones. Stick the electrode pads back onto the protective foil after use and store them in the storage bag to prevent them from drying out. This retains the adhesive force for a longer period.
8. Do not use detergent to clean the electrode pads before and after use to avoid damaging the adhesion of the electrode pads.
9. The electrode pads must always be touched with clean hands, it is recommended to replace the electrode pads if they become dirty.

#### Caution:

- 1) Before applying the electrode, it is recommended for users to wash and degrease the skin, and then dry it.
- 2) Never remove the electrode from the skin while the device is still on.
- 3) Only use the electrode pads provided by the manufacturer. Usage of other companies' products could result in injuries to the user. The electrodes in such way the affected muscles are also enclosed by the electrodes.

#### **6.14 Where do I attach electrode pads?**

1. Each person reacts differently to electric nerve stimulation. Therefore, the placement of the electrodes may deviate from the standard. If application is not successful, contact your physician to find out which placement techniques are best for you.
2. Do not use any adhesive electrodes with a size smaller than those the original manufacturer attached. Otherwise the current density may be too high and cause injuries.
3. The size of the adhesive pads may not be changed, e.g. by clipping off parts of them.
4. Make sure that the region radiating the pain is enclosed by the electrodes. In case of painful muscle groups, attach the electrodes in such a way that the affected muscles are also enclosed by the electrodes.

#### **7. RECOMMENDATIONS FOR USE**

##### **Usage advice for TENS:**

- 1) If you feel the output intensity too strong, you can press [-] button to decrease it;
- 2) If you don't feel any discomfort during the treatment, we advise you to use the device until the session ends. Normally, the pain relief occurs after 5~10 mins treatment;
- 3) Normally, we advise 1~2 treatments per day and one week as a period of treatment;
- 4) After a period of treatment, if the pain relief is not achieved or the pain gets even worse, please consult your doctor.

##### **Usage advice for EMS:**

- 1) Place the electrodes on the body part you want to treat.
- 2) 1~2 treatment per day
- 3) We advise you to use the device for one session per time. If you feel discomfort during treatment, you can either pause the session or decrease the intensity of the output.

##### **Usage advice for MASSAGE**

- 1) Place the electrodes on the body part you want to treat.
- 2) 1~2 treatment per day
- 3) If you feel discomfort during treatment, you can either pause the session or decrease the intensity of the output.
- 4) Massages are mainly used for people whose muscle has stiffness, knotty muscles, tight feeling.

## 8. PROGRAMS AND RELATED INDICATIONS

Below you can find the list of programs and their features. Should the user find a program particularly pleasing and/or satisfactory, he could well also treat other zones with the same programme.

After choosing the treatment, place the electrodes on the area to be treated or take a cue from the illustrations on pages 62 to 65. Or consult tesmed.com under "HELPFUL TIPS" and click on "Positioning electrodes"

Mode	Program	Pulse rate (Hz)	Pulse width (uS)	Treatment time (Min)	Remark	Body part
TENS	1	80-120	120-100	30 (5-90)	Modulation	CERVICAL PAIN
	2	80-100	100	30 (5-90)	Modulation	SHOULDER
	3	2-60	260-160	30 (5-90)	Modulation	ARM
	4	100	150	30 (5-90)	Continue	HAND
	5	10-60	200	30 (5-90)	Modulation	BACK
	6	120	55	30 (5-90)	Continue	ABDOMEN
	7	100	150	30 (5-90)	Burst	WRIST
	8	6-50	200	30 (5-90)	Modulation	LEG
	9	80	180	30 (5-90)	Continue	ANKLE
	10	120	100-120	30 (5-90)	Modulation	(ELBOW/HIP /KNEE/FOOT)
	U1	60 (2-100)	180 (100-300)	30 (5-90)	Continue	ALL BODY PART
EMS	1	30	200	30 (5-90)	Synchronous	BACK OF THE NECK
	2	45	200	30 (5-90)	Synchronous	SHOULDER
	3	50	150	30 (5-90)	Synchronous	ARM
	4	4	200	30 (5-90)	Continue	HAND
	5	60	200	30 (5-90)	Synchronous	BACK
	6	20	200	30 (5-90)	Synchronous	ABDOMEN
	7	30	150	30 (5-90)	Synchronous	WRIST
	8	80	200	30 (5-90)	Synchronous	LEG
	9	20	200	30 (5-90)	Synchronous	ANKLE
		U1	55 (20-100)	200 (100-300)	30 (5-90)	Synchronous

Mode	Program	Pulse rate (Hz)	Pulse width (µs)	Treatment time (Min)	Remark	Body part
M A S S A G E	<b>MESSAGE KNEAD</b> 					
	1	8-24	120-300	30	Modulation	BACK OF THE NECK
	2	16-32	120-280	30	Modulation	SHOULDER
	3	40-56	120-250	30	Modulation	ARM
	4	10-26	120-200	30	Modulation	HAND
	5	45-61	120-300	30	Modulation	BACK
	6	20-36	120-300	30	Modulation	ABDOMEN
	7	6-22	120-280	30	Modulation	WRIST
	8	15-31	120-280	30	Modulation	LEG
	9	4-20	120-150	30	Modulation	ANKLE
	10	22-38	120-200	30	Modulation	HIP
	11	12-28	120-250	30	Modulation	ELBOW
	12	30-46	120-280	30	Modulation	KNEE
	13	40-56	120-300	30	Modulation	FOOT
	<b>MESSAGE RUB</b> 					
	1	8-34	120-300	30	Modulation	BACK OF THE NECK
	2	16-42	120-280	30	Modulation	SHOULDER
	3	40-66	120-250	30	Modulation	ARM
	4	10-36	120-200	30	Modulation	HAND
	5	45-71	120-300	30	Modulation	BACK
	6	20-46	120-300	30	Modulation	ABDOMEN
	7	6-32	120-280	30	Modulation	WRIST
	8	15-41	120-280	30	Modulation	LEG
	9	4-30	120-150	30	Modulation	ANKLE
	10	22-48	120-200	30	Modulation	HIP
	11	12-38	120-250	30	Modulation	ELBOW
12	30-56	120-280	30	Modulation	KNEE	
13	40-66	120-300	30	Modulation	FOOT	

Mode	Program	Pulse rate (Hz)	Pulse width (uS)	Treatment time (Min)	Remark	Body part
M A S S A G E	<b>MESSAGE KNOCK</b> 					
	1	8-30	120-300	30	Modulation	BACK OF THE NECK
	2	16-38	120-280	30	Modulation	SHOULDER
	3	40-62	120-250	30	Modulation	ARM
	4	10-32	120-200	30	Modulation	HAND
	5	45-67	120-300	30	Modulation	BACK
	6	20-42	120-300	30	Modulation	ABDOMEN
	7	6-28	120-280	30	Modulation	WRIST
	8	15-37	120-280	30	Modulation	LEG
	9	4-26	120-150	30	Modulation	ANKLE
	10	22-44	120-200	30	Modulation	HIP
	11	12-34	120-250	30	Modulation	ELBOW
	12	30-52	120-280	30	Modulation	KNEE
	13	40-62	120-300	30	Modulation	FOOT
	<b>MESSAGE TAP</b> 					
	1	8	300	30	Modulation	BACK OF THE NECK
	2	16	280	30	Modulation	SHOULDER
	3	40	250	30	Modulation	ARM
	4	10	200	30	Modulation	HAND
	5	45	300	30	Modulation	BACK
	6	20	300	30	Modulation	ABDOMEN
	7	6	280	30	Modulation	WRIST
	8	15	280	30	Modulation	LEG
	9	4	150	30	Modulation	ANKLE
	10	22	200	30	Modulation	HIP
	11	12	250	30	Modulation	ELBOW
12	30	280	30	Modulation	KNEE	
13	40	300	30	Modulation	FOOT	

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#### **Massage KNEAD:**

It feels like Knead, using series of medium rate pulsing sensations to mimic massage, which likes to manipulate or massage with a kneading motion kneaded muscles.

#### **Massage RUB:**

It feels like Rub, using series of high rate pulsing sensations to mimic hands rubbing, which likes to press against muscles each other with a circular or up-and-down repeated movement.

#### **Massage KNOCK:**

It feels like Knock, using series of medium rate pulsing sensations to mimic massage, which likes to hit muscles with the knuckles of your hand.

#### **Massage TAP:**

It feels like Tap, using series of low rate tapping sensations, which likes to hit muscles lightly.

## 9. CLEANING AND MAINTENANCE

Fully comply with the following necessary daily maintenance requirements to make sure the device is intact and guarantee its long-term performance and safety.

### 9.1 Cleaning and care for the device

1. Pull the electrodes out of the stimulator, clean the device with a soft, slightly damp cloth. In case of heavier dirt build-up, you may also apply a mild detergent.
2. Do not expose the TE-880R PLUS stimulator to moisture or dampness. And do not hold the TE-880R PLUS stimulator under running water, nor submerge it in water or other liquids.
3. The TE-880R PLUS stimulator is sensitive to heat and may not be exposed to direct sunlight. And do not place it on hot surfaces.
4. Clean the surface of the electrode pads carefully with a damp cloth. Make sure the device is turn off.
5. For reasons of hygiene, each user should use his/her own set of electrodes.
6. Do not use any chemical cleaners or abrasive agents for cleaning.
7. If water has penetrated the device, do not use it and contact the service department by e-mailing [info@tesmed.com](mailto:info@tesmed.com)
8. Do not clean the device during treatment. Be sure that the device is turned off before cleaning.

### 9.2 Maintenance

1. The manufacturer didn't authorize any maintenance agencies abroad. If your device has any problems, please contact [info@tesmed.com](mailto:info@tesmed.com). The manufacturer will not be responsible for the results of maintenance or repairs by unauthorized persons.
2. The user must not attempt any repairs to the device or any of its accessories. Please contact [info@tesmed.com](mailto:info@tesmed.com) for repair.
3. Opening of the equipment by unauthorized agencies is not allowed and will terminate any claim to warranty.

Each product in manufacturing has been inspected through the systematic validation. The performance is stable and does not need to undertake calibration and validation.

## 10. TROUBLESHOOTING

Malfunction	Common reasons	Countermeasure
No display	<ol style="list-style-type: none"><li>1. The battery is exhausted</li></ol>	<ol style="list-style-type: none"><li>1. Charge in time</li></ol>
No sensation of stimulation or weak stimulation	<ol style="list-style-type: none"><li>1. The electrode does not connect well to the skin.</li><li>2. If the connection between electrode connects well to the stimulator.</li><li>3. The battery is used up.</li><li>4. The skin is too dry.</li></ol>	<ol style="list-style-type: none"><li>1. Check the adhesion of the electrode if necessary replace it</li><li>2. Check the connection.</li><li>3. Charge.</li><li>4. Wipe the electrode and the skin with a wet cotton cloth.</li></ol>
Automatic halt in the treatment	<ol style="list-style-type: none"><li>1. The electrode loses connection with the skin.</li><li>2. The battery is used up.</li></ol>	<ol style="list-style-type: none"><li>1. Check the adhesion of the electrode if necessary replace it</li><li>2. Charge</li></ol>

Malfunction	Common reasons	Countermeasure
Rash or tickle on the skin occurs in the treatment	<ol style="list-style-type: none"><li>1. The treatment time lasts too long.</li><li>2. The electrode does not stick well to the skin.</li><li>3. The interface of the electrodes is dirty or dry.</li><li>4. The skin is sensitive to the electrode.</li></ol>	<ol style="list-style-type: none"><li>1. Do the treatment once a day and shorten the treatment time.</li><li>2. Check and stick the electrode well.</li><li>3. Try cleaning the electrode before use and replace it if necessary</li><li>4. Check your allergic history. Please change the sticking place or shorten the treatment time. If your skin is over-sensitive, you should stop the treatment or go to see a doctor.</li></ol>

## **11. STORAGE**

### **11.1 Storing the Electrode Pads and Lead Wires**

1. Turn the device off and remove the lead wires from the unit.
2. Remove the electrodes from your body and disconnect the lead wires from the electrodes.
3. Place the electrodes onto the plastic film and the store into the sealed package.
4. Wrap the lead wires and store into the sealed package.

### **11.2 Storing the Unit**

1. Place the unit, electrodes, lead wires and manual back into the gift box. Store the box in a cool, dry place,  $-10^{\circ}\text{C} \sim 55^{\circ}\text{C}$ ; 10% ~ 90% relative humidity.
2. Do not keep in places that can be easily reached by children

## 12. DISPOSAL

Spent batteries do not belong to the household wastes. Disposal of the battery according to the current regulations. As a consumer, you have the obligation to dispose of batteries correctly.

Consult your municipal authority or your dealer for information about disposal.

At the end of the product lifecycle, do not throw this product into the normal household garbage, but bring it to a collection point for the recycling of electronic equipment. Obsolete electrical and electronic equipment may have potentially harmful effects on the environment. Incorrect disposal can cause toxins to build up in the air, water and soil and jeopardize human health.

### 13. ELECTROMAGNETIC COMPATIBILITY (EMC) TABLES

Guidance and manufacture's declaration - electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
R F emissions CISPR11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
R F emissions CISPR11	Class B	The device is suitable for use in all establishments including those directly connected to the public low-voltage power supply network that supplies to buildings power used for domestic purposes
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Not applicable	
Guidance and manufacture's declaration — electromagnetic immunity		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environment.		

Testmed copyright - It is absolutely forbidden the partial or total reproduction of any present content

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8kV direct & indirect contact; ±15kV air discharge	±8kV direct & indirect contact; ±15kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	Not applicable	not applicable (for INTERNALLY POWERED ME EQUIPMENT)
Surge IEC 61000-4-5	± 1 kV line(s) to line(s)	Not applicable	not applicable (for INTERNALLY POWERED ME EQUIPMENT)
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles  70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 sec	Not applicable	not applicable (For INTERNALLY POWERED ME EQUIPMENT)

Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	10V/m	10V/m	Power frequency magnetic fields should be at levels characteristic of a typical location in typical commercial or hospital environment.
---	-------	-------	---

NOTE UT is the a.c. mains voltage prior to application of the test level.

#### Guidance and manufacture's declaration – electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Radiated RF IEC	10V/m & table 9	10V/m & table 9	<p>Portable and mobile RF communications equipment should be used not closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.167 \sqrt{P}$ <p>to 80 MHz</p> $d = 2.333 \sqrt{P}$ <p>to 800 MHz</p> <p>a 2,5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended</p>

Tested copyright - It is absolutely forbidden the partial or total reproduction of any present content.

		<p>separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. b Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
--	--	--

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people

a. Field strengths from fixed transmitters, such as the base station for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcasts and TV broadcasts cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey must be considered. If the field strength measured at the location where the device is used exceeds the applicable RF compliance level above, the device must be observed to verify normal operation.

If abnormal performance is observed, additional measures, such as reorienting or relocating the device, may be necessary.

b. In the frequency range between 150 kHz and 80 MHz, field strengths must be less than  $[V_i] \text{ V/m}$ .

**Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment (Table 9)**

Test frequency (MHz)	Band <sup>a)</sup> (MHz)	Service <sup>a)</sup>	Modulation <sup>a)</sup>	Maximum power (W)	Distance (m)	Immunity Test Level (V/m)
385	380-390	TETRA 400	Pulse modulation <sup>b)</sup> 18Hz	1.8	0.3	27
450	430-470	GMRS 460, FRS 460	FM <sup>c)</sup> ±5kHz deviation 1kHz sine	2	0.3	28
710 745 780	704-787	LTE Band 13, 17	Pulse modulation <sup>b)</sup> 217Hz	0.2	0.3	9
810 870 930	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation <sup>b)</sup> 18Hz	2	0.3	28
1720 1845 1970	1700-1990	GSM1800; CDMA 1900; GSM 1900; DECT; LTE Band 1,3,4,25; UMTS	Pulse modulation <sup>b)</sup> 217Hz	2	0.3	28
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation <sup>b)</sup> 217Hz	2	0.3	28
5240 5500 5785	5100-5800	WLAN 802.11 a/n	Pulse modulation <sup>b)</sup> 217Hz	0.2	0.3	9

Tested copyright - it is absolutely forbidden the partial or total reproduction of any present content

NOTE If it is necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

- a) For some services, only the uplink frequencies are included.
- b) The carrier shall be modulated using a 50 % duty cycle square wave signal.
- c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because it does not represent actual modulation. It would be the worst case.

## 14. NORMALIZED SYMBOLS

	WEEE Symbol		Applied part of type BF
	Refer to instruction manual	<b>IP22</b>	IP classification
<b>LOT</b>	Batch code	<b>EC REP</b>	Authorized representative in the European Community
	Manufacturer		Date of manufacture
	Fragile, handle with care		Keep away from sunlight
	Keep dry		Temperature limit
	Use-by date	<b>SN</b>	Serial number
	Humidity limitation		Atmospheric pressure limitation
	Caution		Recycling symbol
	Users of the artificial pacemaker are prohibited from using the device	<b>CE</b> 2460	CE mark
	Packaging material cycle mark	<b>MD</b>	Medical device
<b>UDI</b>	Unique device identifier		Importer

## 15. WARRANTY

### HOW THE WARRANTY IS EXERCISED.

To exercise the rights of the legal guarantee, the consumer can contact, with a receipt proving the purchase made, at the point of sale where he purchased the product. The store where the purchase was made cannot refuse to take delivery of the non-compliant product for repair or replacement.

If the purchase was made through internet sites, the user must keep the shipping confirmation received by email and write directly to the seller on the website on which the product was purchased to request repair or replacement under warranty.

For online purchases made through the website [www.tesmed.it](http://www.tesmed.it) as well as [www.tesmed.com](http://www.tesmed.com), it is sufficient to keep the shipping confirmation email and / or the order confirmation email.

For the conditions of the guarantee, visit the [tesmed.com](http://tesmed.com) website in the section dedicated to legal notes, conditions of sale and the guarantee.

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© 2005 Tesmed - All rights reserved. This is a reproduction for sale in Italy of a quality controlled product.
© 2004 by Elsevier GmbH. All rights reserved. No part of this publication may be reproduced, stored in a retrieval system, or transmitted, in any form or by any means, electronic, mechanical, photocopying, recording, or otherwise, without the prior written permission of Elsevier GmbH.
© 2008 Tesmed - all rights reserved. All trademarks are the property of their respective owners.
Shenzhen Roundwhale Technology Co.,Ltd.
Shanghai International Holding Corp. GmbH(Europe)