PROFESSIONAL CARE

designed by **BRAUN**

OXYJET

Deutsch	6
English	11
Français	15
Español	20
Português	24
Italiano	28
Nederlands	32
Dansk	36
Norsk	40
Svenska	44
Suomi	48

DE AT	00 800 27 28 64 63
CH	08 44 - 88 40 10
GB	0 800 731 1792
IE	1 800 509 448
FR	0 800 944 803 (service et appel gratuits)
BE	0 800 14 592
LU	800 21172
ES	900 814 208
PT	808 20 00 33
IT	800 124 600
NL	0800-4 45 53 88
DK	70 15 00 13
NO	22 63 00 93
SE	020 - 21 33 21
FI	020 377 877

Internet:

www.oralb.com
www.service.braun.com

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English

Welcome to Oral-B!

Before operating this appliance, please read these instructions and save this manual for future reference.

IMPORTANT

- The installation instructions shall state that parts that have to be fixed must be fixed so that they cannot fall into water, unless they are of IPX7 construction.
- Periodically check the entire product/cord for damage. A damaged or non-functioning unit should no longer be used. Do not modify or repair the product. This may cause fire, electric shock or injury.
- If the product/ cord is damaged, take it to an Oral-B Service Centre.
- Never insert any object into any opening of the appliance.
- Oral irrigators can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge under supervision.
- Oral irrigators shall be used only with liquids defined in the manufacturer's instructions.
- Do not use iodine, saline (salt solution), or water insoluble concentrated essential oils in the appliance. Use of these can reduce product performance and will shorten the life of the product. Use of these products may void warranty.
- Use this product only for its intended use as described in this manual. Do not use attachments which are not recommended by the manufacturer.
- The appliance must be mounted in a way that it cannot fall into water or drop to the floor. The plug should always be accessible.
- Only use charger provided with your appliance.
- Toothbrushes can be used by children and persons with reduced physical,

sensory or mental capabilities or lack of experience and knowledge, if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children shall not play with the appliance. Cleaning and user maintenance shall not be made by children without supervision.

WARNING

- Do not place the charger, smart plug or charging travel case in water or liquid or store where it can fall or be pulled into a tub or sink. Do not reach for it when fallen into water. Unplug immediately.
- Avoid sharp bending, jamming or squeezing of the cord.
- This appliance contains batteries that are non-replaceable. Do not open and disassemble device. For battery recycling please dispose whole device according to local environmental regulations. Opening the handle will destroy the appliance and invalidate the warranty.
- For tooth brushing only.
- When unplugging, always pull the power plug instead of the cord. Do not touch the power plug with wet hands. This can cause an electric shock.
- Do not chew or bite on brush. Do not use as a toy or without adult supervision.
- Usage by children under age 3 is not recommended.
- This toothbrush is a personal care device and is not intended for use on multiple patients in a dental practice or institution.
- For hygienic reasons we recommend the use of individual brush handles and brush heads for each family member.
- Do not apply stickers on handle.
- Small parts may come off, keep out of reach of children.
- This irrigator is a personal care device and is not intended for use on multiple patients in a dental practice or institution.
- Do not chew or bite on nozzle. Do not use as a toy or without adult supervision.
- For water flossing only.
- Do not direct water under the tongue, into the ear, nose or other delicate areas. This product is capable of producing pressures that may cause serious damage in these areas.
- Do not direct water into the nose or ear. The potentially deadly amoeba, *Naegleria fowleri*, may be present in some tap water or unchlorinated well water and may be fatal if directed into these areas.

To avoid brush head breakage which may generate a small part choking hazard or damage teeth:

- Prior to each use, ensure the brush head fits well. Discontinue use of the toothbrush, if brush head no longer fit properly. Never use without a brush head.
- If the toothbrush handle is dropped, the brush head should be replaced before the next use even if no damage is visible.
- Clean brush head properly after each use (see «Cleaning recommendations» paragraph). Proper cleaning ensures the safe use and functional life of the toothbrush.
- Replace brush (head) every 3 months or sooner if brush head becomes worn.

Description

- a Nozzle
- b Nozzle switch
- c Nozzle compartment
- d Nozzle release button
- e On/off switch
- f Water intensity control dial
- g Handle switch (water stop and go)
- h Irrigator handle
- i Water container (600 ml)
- j Charger
- k Charge/ low-charge indicator light
- l Toothbrush handle
- m Toothbrush on/off button (mode selection)
- n Brush head
- o Cord storage
- p Wall holder

Specifications

For voltage specifications refer to the bottom of the charging unit.

Noise level toothbrush: ≤ 65 dB (A)

Noise level Oral irrigator: 76 dB (A)

A - Using your Oral irrigator

Your Oral irrigator comes with a convenient wall holder. It can be mounted easily using the attached items (picture 1).

Before using for the first time, fill the water container (i) and allow the water to run through completely by holding the irrigator handle (h) over a washbasin.

1. Attach one of the irrigator nozzles (a) to the handle until it locks into place (picture 1). Fill the water container with luke warm water.
2. Turn on the irrigator by pressing the on/off switch (e) (picture 3). (When doing so, the handle switch should be in the lower position = water stop).
3. Select the water pressure with the pressure control dial (f) (1 = gentle, 5 = strong) (picture 4). We recommend that you start with setting 1, especially if your gums tend to bleed or if children are using the Oral irrigator.
4. Bend over a washbasin and guide the nozzle to your teeth and gums. Slide the handle switch (g) up to start the water stream (picture 5). Keep your mouth slightly open to allow the water to run into the washbasin.

12

Slide the handle switch up (g) for a continuous water stream.

Slide the handle switch down and press into it, to use the on-demand function (picture 6).

Caution: Turn off the Oral irrigator by switching off the on/off switch once the water container is empty. Longer operation of the Oral irrigator without water may damage the appliance.

OxyJet nozzle

Nozzle functions:

You can choose between two different nozzle functions by sliding the nozzle switch (b) down or up. Before doing so, turn off the water stream by sliding the handle switch down:

Rotating stream for general cleaning and gum massage (nozzle switch in lower position): The rotation of the water stream and its micro air-bubbles are designed for easy penetration into the interdental space and gentle gum massage (picture 8).

Focused stream type (nozzle switch in upper position):

The focused water stream with its micro-bubbles are designed for targeted cleaning (picture 9).

WaterJet nozzle

Jet functions:

You can choose between different stream types by click-in and click-out the blue center of the nozzle. Before doing so, turn off the water stream by sliding the handle switch down:

When the blue center is out (default position) you will get a high pressure single precision stream (picture 10).

As soon as you click the blue center in, you will get a multi-jet stream, designed for great reach and massage effect of the gums and the gumline (picture 11).

Maximal operation time: 15 minutes

Cooling off time: 2 hours

After use

Always empty the water container entirely. Keep the Oral irrigator switched on until the water stops running through the nozzle in order to prevent accumulation of bacteria. As soon as no more water flows out of the nozzle, switch the device off (picture 3). Press the nozzle release button (d) (picture 7) to take off the nozzle from the Oral irrigator. Always wipe the handle and the Oral Irrigator dry; place the handle back onto its base and store the Oral irrigator in the nozzle compartment (c).

Cleaning recommendation

After use, rinse the nozzle under running water. Regularly clean the housing with a damp cloth. The water container and the nozzle compartment are removable and dishwasher-safe.

B - Using your toothbrush

Charging and Operating

This toothbrush has a waterproof handle, is electrically safe and designed for use in the bathroom.

- Plug the charging unit (e) into an electrical outlet and place the toothbrush handle (c) on the charging unit (picture 1).
- The green charge indicator light (d) flashes while the handle is being charged. Once it is fully charged, the light turns off. A full charge takes about 22 hours. Note: After deep discharge, the green light might not flash immediately; it can take up to 10–15 minutes.
- If the rechargeable battery is low the red charge indicator light (d) is flashing for some seconds when turning your toothbrush on/off.
- For everyday use, the handle can be stored on the plugged-in charging unit to maintain it at full power. Overcharging is impossible.
- To maintain the maximum capacity of the rechargeable battery, unplug the charging unit and fully discharge the handle with regular use at least every 6 months.

Using the toothbrush

Brushing technique

- Wet brush head and apply any kind of toothpaste. To avoid splashing, guide the brush head to your teeth before switching on the appliance (picture 2).
- When brushing your teeth with one of the Oral-B oscillating-rotating brush heads guide the brush head slowly from tooth to tooth, spending a few seconds on each tooth surface.
- With any brush head start brushing the outsides, then the insides and finally the chewing surfaces (picture 4). Brush all four quadrants of your mouth equally (picture 4). Do not press too hard, simply let the brush do all the work. You may also consult your dentist or dental hygienist about the right technique for you.

In the first few days of using any electric toothbrush, your gums may bleed slightly. In general, bleeding should stop after a few days. Should it persist after 2 weeks, please consult your dentist or dental hygienist.

If you have sensitive teeth and/or gums, Oral-B recommends using the «Sensitive» mode (optional in combination with an Oral-B «Sensitive Clean» brush head).

Brushing modes (depending on model)

Your toothbrush offers different brushing modes:

«**Daily Clean**» – Standard mode for daily cleaning

«**Sensitive**» – Gentle, yet thorough cleaning for sensitive areas

«**Whitening**» – Polishing for occasional or everyday use

How to switch between modes

Your toothbrush automatically starts in the «Daily Clean» or «Sensitive» mode (depending on model). To switch to other modes, successively press the on/off button. Switch off the toothbrush by holding the on/off button until the motor stops.

Pressure Sensor

To promote optimal brushing your toothbrush has a pressure control feature installed. If too much pressure is applied, the movement of the brush head will continue but its pulsation will stop. In addition, you will also hear a different sound while brushing.

«Professional» timer:

A short stuttering sound at 30-second intervals reminds you to brush all four quadrants of your mouth equally (picture 3). A long stuttering sound indicates the end of the professionally recommended 2-minute brushing time. If the toothbrush is turned off during brushing, the reached brushing time will be memorised for 30 seconds. When pausing longer than 30 seconds the timer resets.

Brush heads

Oral-B offers you a variety of different Oral-B brush heads that fit your Oral-B toothbrush handle.

Oral-B CrossAction brush head

Our most advanced brush head. Angled bristles for a precise clean. Lifts and powers away plaque.

Oral-B Sensitive Clean brush head

features a combination of regular bristles to clean teeth surfaces and ultrathin bristles for a gentle brushing experience on the gumline.

Most Oral-B brush heads feature INDICATOR® bristles to help you monitor brush head replacement need. With thorough brushing, twice a day for two minutes, the colour will fade approximately within 3 months, indicating the need to exchange your brush head. If the bristles splay before the colour recedes, you may exert too much pressure on teeth and gums (picture 4).

We do not recommend using the Oral-B «FlossAction» or Oral-B «3D White» brush head with braces. You may use the Oral-B «Ortho» brush head, specifically designed to clean around braces and wires.

Cleaning recommendations

After use, rinse the brush head for several seconds under running water with the handle switched on. Switch off and remove the brush head. Rinse both parts separately under running water and wipe them dry before reassembling the toothbrush.

Unplug the charger before cleaning. Brush head holder (f) is dishwasher-safe. Charger (e) and travel case (g) should be cleaned with a damp cloth only (picture 5). Store the clean and dry toothbrush/brush heads in the travel case.

Subject to change without notice.

Environmental Notice

Product contains batteries and/or recyclable electric waste. For environment protection do not dispose in household waste, but for recycling take to electric waste collection points provided in your country.

- Oral-B cannot predict the long-term effect of non-Oral-B replacement brush heads/nozzles on handle/irrigator wear.

All Oral-B replacement brush heads/nozzles carry the Oral-B logo and meet Oral-B's high quality standards. Oral-B does not sell replacement brush heads/nozzles or handle/irrigator parts under any other brand name.

Guarantee

We grant 2 years guarantee on the product commencing on the date of purchase.

Within the guarantee period we will eliminate, free of charge, any defects in the appliance resulting from faults in materials or workmanship, either by repairing or replacing the complete appliance as we may choose. Depending on availability replacement of a unit might result in a different colour or an equivalent model.

This guarantee extends to every country where this appliance is supplied by Braun / Oral-B or its appointed distributor. This guarantee does not cover: damage due to improper use, normal wear or use, especially regarding the brush heads/nozzles, as well as defects that have a negligible effect on the value or operation of the appliance.

The guarantee becomes void if repairs are undertaken by unauthorized persons and if original Braun / Oral-B parts are not used.

To obtain service within the guarantee period, hand in or send the complete appliance with your sales receipt to an authorized Oral-B Braun Customer Service Centre.

This guarantee in no way affects your rights under legal/ statutory law.

Warranty replacement brush heads/nozzles

Oral-B warranty will be voided if electric rechargeable toothbrush handle/irrigator damage is found to be attributed to the use of non-Oral-B replacement brush heads/nozzles.

Oral-B does not recommend the use of non-Oral-B replacement brush heads/nozzles.

- Oral-B has no control over the quality of non-Oral-B replacement brush heads/nozzles. Therefore, we cannot ensure the cleaning performance of non-Oral-B replacement brush heads/nozzles, as communicated with the electric rechargeable handle/irrigator at time of initial purchase.
- Oral-B cannot ensure a good fit of non-Oral-B replacement brush heads/nozzles.

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