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directly, e.g. 32.1 °C (89,9°F). If temperatures over 37.8°C (100°F) are measured, a fever alarm (30 short beeps in 10 seconds) sounds. If temperatures over 43.0 °C (109.4°F) are measured, the symbol for »Hi °C (Hi °F) « appears. The measurement time is very quick compared to conventional thermometers and is only 10 seconds for the rectal measurement. More accurate measurements are possible by extending the measurement time beyond the acoustic signal. For measurements in the armpits the recommended measurement time should be maintained in spite of the acoustic signal. To prolong the life of the battery, please switch off the thermometer by briefly pressing the On/ Off key after the measurement. Otherwise the thermometer will automatically switch off after 10 minutes.

Memory
The last measured value will be stored automatically. This value is displayed when the On/Off button is pressed. The memory symbol is also shown beside the measured value, e.g. 37.5 °CM (99.5 °FM).

Changing the measuring unit
The temperature can be displayed in degrees Celsius (°C) and degrees Fahrenheit (°F). To set the measuring unit press and hold the On/Off button for longer than 2 seconds when switching on until the measuring unit changes. The saved value is displayed in the new measuring unit.

Changing the battery
The clinical thermometer contains a long-life alkali-manganese battery (LR 41, 1.55 V), which will last for approx. 2 years. The battery has to be replaced when "T" appears in the bottom right corner of the display. Press down the battery cover at the end of the device with your fingers and then pull it firmly back. Remove the old battery and replace it with a battery of the same type. The + sign on the battery must point upwards. Then put the battery cover back onto the case. Used batteries do not belong in the regular domestic refuse. Dispose of them through your local dealer for electric goods or your local depot for special waste. You are legally obliged to dispose of the batteries.

Note: Batteries containing pollutant substances are marked as follows: Pb = Battery contains lead, Cd = Battery contains cadmium, Hg = Battery contains mercury.

Cleaning/Disinfection
The thermometer is 100% water-proof. For the purpose of cleaning or disinfection, it can be dipped into water or disinfectant without any difficulties.
Under no circumstances should the thermometer be boiled! Use one

directly, e.g. 32.1 °C (89,9°F). If temperatures over 37.8°C (100°F) are measured, a fever alarm (30 short beeps in 10 seconds) sounds. If temperatures over 43.0 °C (109.4°F) are measured, the symbol for »Hi °C (Hi °F) « appears. The measurement time is very quick compared to conventional thermometers and is only 10 seconds for the rectal measurement. More accurate measurements are possible by extending the measurement time beyond the acoustic signal. For measurements in the armpits the recommended measurement time should be maintained in spite of the acoustic signal. To prolong the life of the battery, please switch off the thermometer by briefly pressing the On/ Off key after the measurement. Otherwise the thermometer will automatically switch off after 10 minutes.

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of the following disinfectants in accordance with the manufacturer's instructions. Other cleaning agents or methods could impair the function of the thermometer or damage it.

Name of the agent: Isopropanol 70% Pharmacy; Gigasept FF; Lysoformin; Aseptisol.

Storage/Disposal
When you are not using the thermometer, please retain it in the original protective cover. The device must be tested for accuracy (calibrated) after 2 years. Please send the complete device to the dealer's or service address.

Please dispose of the device in accordance with the directive 2002/96/EG – WEEE (Waste Electrical and Electronic Equipment). If you have any queries, please refer to the local authorities responsible for waste disposal.

Technical data
Type: Maximum thermometer type FT 15/1
Range of measurement: 32 °C (89,6°F) to 42.9 °C (109,2°F)
Accuracy of of measurement: ±0.1 °C in a water bath with a temperature of 35,5 °C (95,9°F) to 42,0 °C (107,6°F)
±0.2 °C in a water bath with a temperature of 32,0 °C (89,6°F) to 35,4 °C (95,7°F)
±0.2 °C in a water bath with a temperature of 42,1 °C (107,8°F) to 42,9 °C (109,2°F)

Ambient temperature: 10 °C (50°F) to 40 °C (122°F), with 30% to 85% relative humidity

Storage temperature: -10 °C (14°F) to +60 °C (140°F), with 25% to 90% relative humidity

Reference of standardisation: EN 12470-3: 2000 Clinical Thermometers

Explanation of symbols: Attention – read the operating instructions.

Application part of Type BF

of the following disinfectants in accordance with the manufacturer's instructions. Other cleaning agents or methods could impair the function of the thermometer or damage it.

Name of the agent: Isopropanol 70% Pharmacy; Gigasept FF; Lysoformin; Aseptisol.

Storage/Disposal
When you are not using the thermometer, please retain it in the original protective cover. The device must be tested for accuracy (calibrated) after 2 years. Please send the complete device to the dealer's or service address.

Please dispose of the device in accordance with the directive 2002/96/EG – WEEE (Waste Electrical and Electronic Equipment). If you have any queries, please refer to the local authorities responsible for waste disposal.

Technical data
Type: Maximum thermometer type FT 15/1
Range of measurement: 32 °C (89,6°F) to 42.9 °C (109,2°F)
Accuracy of of measurement: ±0.1 °C in a water bath with a temperature of 35,5 °C (95,9°F) to 42,0 °C (107,6°F)
±0.2 °C in a water bath with a temperature of 32,0 °C (89,6°F) to 35,4 °C (95,7°F)
±0.2 °C in a water bath with a temperature of 42,1 °C (107,8°F) to 42,9 °C (109,2°F)

Ambient temperature: 10 °C (50°F) to 40 °C (122°F), with 30% to 85% relative humidity

Storage temperature: -10 °C (14°F) to +60 °C (140°F) for a humidity relative of 25% to 90%

Reference of standardisation: EN 12470-3 : 2000 Clinical Thermometers

Explanation of symbols: Attention – read the operating instructions.

Application part of Type BF

of the following disinfectants in accordance with the manufacturer's instructions. Other cleaning agents or methods could impair the function of the thermometer or damage it.

Name of the agent: Isopropanol 70% Pharmacy; Gigasept FF; Lysoformin; Aseptisol.

Storage/Disposal
When you are not using the thermometer, please retain it in the original protective cover. The device must be tested for accuracy (calibrated) after 2 years. Please send the complete device to the dealer's or service address.

Please dispose of the device in accordance with the directive 2002/96/EG – WEEE (Waste Electrical and Electronic Equipment). If you have any queries, please refer to the local authorities responsible for waste disposal.

Technical data
Type: Maximum thermometer type FT 15/1
Range of measurement: 32 °C (89,6°F) to 42.9 °C (109,2°F)
Accuracy of of measurement: ±0.1 °C in a water bath with a temperature of 35,5 °C (95,9°F) to 42,0 °C (107,6°F)
±0.2 °C in a water bath with a temperature of 32,0 °C (89,6°F) to 35,4 °C (95,7°F)
±0.2 °C in a water bath with a temperature of 42,1 °C (107,8°F) and 42,9 °C (109,2°F)

Ambient temperature: 10 °C (50°F) and 40 °C (122°F) for a humidity relative of 30% to 85%

Storage temperature: -10 °C (14°F) to +60 °C (140°F) for a humidity relative of 25% to 90%

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