currentbody.com

CurrentBody.com Ltd, Q17, Crossley Road, Stockport, Greater Manchester,
SK4 5BB, United Kingdom.

CurrentBody BV, Naritaweg 2, 1437EL, Rozenburg, Netherlands

Manual Version: V1.0
Issue Date: 2022-03-08

CURRENTBODY | skin™

LED LIP PERFECTOR

USER GUIDE

Includes important safety information.
Read all instructions before using device.

Welcome to CurrentBody Skin™ LED Lip Perfector

Established in 2009, CurrentBody is the world's only online retailer dedicated to home-use beauty devices. As the first company to recognise devices as one of the fastest growing trends in the world of beauty, we have unrivalled expertise in selecting at-home devices that are science-backed and clinically proven. We offer our customers the ultimate choice of devices that treat a wide range of concerns, including anti-ageing, cleansing and hair removal.

Underpinning our range of innovative beauty technology is expert education to help our customers choose the right device for them and excellent customer service.

Our highly trained beauty device experts are on hand to offer impartial advice and answer your questions in person, as well as how-to videos and top tips available throughout our website and social media channels.

CurrentBody Skin LED Lip Perfector is a light therapy device for treatment of lip creases and wrinkles, and lip plumping.

WHAT IS IN THE BOX

Lip Care Device

Removable Mouth Tray

User Guide

USB charging cord

Storage bag

CONTRAINDICATIONS

Do not use this device if your skin is sensitive to light or you are currently using medication that may cause skin sensitivity to light. Doing so may result in a photosensitive reaction.

If you are unsure whether your medication causes sensitivity to light, confirm with your physician prior to using this device.

Discontinue use immediately if you experience any sensitivity to using this device.

WARNINGS

1. Do not operate or charge this device if it has come into contact with water.
2. Stop using device immediately if you think it is defective.
3. Do not attempt to take the device apart. Any modification or attempted service will void the warranty.
4. Do not store the device in sun or a location with high temperatures as this may cause damage to plastic casing.
5. Discontinue use and contact Customer Service if device overheats or becomes uncomfortably hot to your skin.
6. Always unplug the device from the outlet after charging is complete.
7. Use only the supplied Universal USB charging cord. Use of an incorrect power cord could damage the device.
8. Do not use if your skin is red or swollen from a previous treatment.
9. Do not use the device to treat skin that is painful to touch, raised (higher than surrounding skin), crusting, or bleeding. This may be a sign of a more serious condition and should be checked by your doctor.
10. This device contains a Lithium Ion Battery. Before disposal, please review your local laws and requirements surrounding Lithium Ion Battery disposal. The preferred method of disposal is through recycling.

IMPORTANT SAFETY INFORMATION

This device is designed for use on adults over 18 years of age and not recommended for use on infants or children.

Do not shine the device lights directly into the eyes. Doing so may cause temporary vision discomfort.

Contact your doctor if you have any medical questions regarding use of this device.

Consult your doctor prior to using this device if you:

Are pregnant or suspect you may be pregnant. Have suspicious cancerous lesions.

Customer Care:

Please email sales@currentbody.com if you have any questions or you need any advice. Do not use or charge the device with a power adapter that is damaged or faulty. Do not harshly bend or abuse the USB charging cord.

DEVICE OVERVIEW

CurrentBody Skin™ Lip Care Device is a light emitting diode (LED) device that harnesses energy to plump lips and treat wrinkles around the lip area. The device incorporates natural light energy in Amber, Red, Deep Red, and Infrared spectrums that are proven to reduce wrinkles, increase circulation, and plump lips.

LIGHT-SENSITIVITY TEST

Test yourself for light-sensitivity prior to initial treatment:

1. Power on your device.
2. Place device on the skin of your forearm for three minutes.
3. After 3 minutes, check the area. If area has turned red, and redness lasts for more than 2 hours, your skin is light-sensitive.
4. If your skin is light-sensitive, return the device to CurrentBody within 14 days for a full refund.

HOW TO USE

Set up:

1. Unpack your CURRENTBODY Skin™ Led Lip Perfector. There are 2 ways to charge your device.

A. From the charging base:

- Place the Led Lip Perfector on the base, ensuring on either side are in contact with the base.
- Plug the USB charging cord into the charging base.
- Plug the USB charging cord into a Universal USB Power Adapter (not included) and then into a mains outlet.

B. By plugging the cord directly into the device and then into a Universal USB Power Adapter (not included). Plug the adapter into a mains outlet.

2. The power indicator light will flash while charging. Charge your device for 1.5 - 2 hours. The device is fully charged when the power indicator light is solid (not flashing). The USB port can be found under the mouthguard. Simply remove and plug in.
3. Place the device in the protective (no gap) storage bag when not in use.
4. Ensure you unplug the device from the USB charger before using. For your safety, this device will not work when it is charging.
5. When the charge is low the indicator light will flash. To recharge the device, follow the instructions above.

Treatment:

1. To use the device, ensure the silicone mouth tray is clean and attach to the device.
2. Place the tray in your mouth so the device is comfortably against your lips.
3. Switch ON. The device will automatically switch off after 3 minutes.
4. You can switch the device off at any time by pressing the ON / OFF button
5. For best results, use daily.

MAINTENANCE AND CLEANING

Maintenance:

1. This device has no serviceable parts. The battery is not accessible and not serviceable. Users should not attempt to access or service the battery as this will damage the device and void the warranty.
2. Device includes only an ON/OFF button. All settings are set to the factory setting.
3. Charge the device only with the supplied USB power charging cord and a universal USB power adapter (not included).
4. The power indicator light will begin to flash when battery levels are low. This indicates the device should be recharged.
5. Compatible Power Supply: 100V/240V.

Cleaning:

1. Turn device off and unplug from any power source.
2. Separate the LED device from the removable mouth tray.
3. Wash mouth tray in lukewarm water and dry.
4. Wipe the device surface with a clean, damp cloth.
5. Allow device to dry completely prior to use.

TROUBLESHOOTING

1. System not turning on:

- Is the device charged? Please try recharging.
- If the device still does not turn on it is faulty. Please contact CurrentBody Customer Service.

2. System shuts off abruptly prior to completion of 3-minute treatment:

- System is designed to shut off automatically after 3 minutes.
- System may shut off automatically after multiple repeated uses if the device exceeds desired temperature. Allow the device to cool for several minutes and try again.
- If the device continues to shut off prior to 3 minutes, please contact CurrentBody customer service as the device is faulty.

3. Some lights appear turned off:

- The Infrared (IR) light spectrum is not visible to the human eye. These lights will look as if they are not on, but they are operating properly. If the red lights are visible, then the device is operating as designed

WARRANTY

Under normal use this product is guaranteed for 2 years after its original purchase date against manufacturing defects in materials or workmanship. This warranty extends only to the original purchaser. If the product should become faulty within the warranty period, contact CurrentBody.

- CurrentBody reserves the right to repair a defective product.
- CurrentBody reserves the right to replace a defective product with the most compatible product currently available.
- This warranty does not cover damage caused by the following:
 - Accident, misuse, abuse or alteration
 - Servicing by unauthorised persons
 - Connection to incorrect power sources
 - Use with unauthorised accessories

CurrentBody shall not be responsible for any incidental, special or consequential damages resulting from the use of this product. All implied warranties, including but not limited to implied warranties of fitness and merchantability are limited in duration to two years from date of original purchase. This warranty gives you specific legal rights, and you may also have other rights that vary from country to country. Some countries do not allow the exclusion or limitation of incidental or consequential damages, or limitations on how long an implied warranty lasts, so some of the above limitations and exclusions may not apply to you. This warranty is given in lieu of all other warranties, written or oral, whether expressed by affirmation, promise, description, drawing, model or sample. Any and all warranties other than this one, whether expressed or implied, including implied warranties of merchantability and fitness for a particular purpose are hereby disclaimed. If the product is defective within the warranty period please contact CurrentBody.

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