Z I I P

HALO

USER GUIDE

CONTENTS

What is in the box?	5
Important safety information	6
Warnings	6
Cautions	8
Contraindications	9
Device Overview	10
How to use	12
Warranty	17
Specification	19
Symbols Explained	20
Safety and EMC	21
Manufacturer information	22
Français	26
Deutsch	44
Español	62
Italiano	80
Svenska	98
Suomi	114
Dansk	130
Norsk	146
繁體中文	162
简体中文	176
日本語	190
한국어	204
Melayu	218
Türkçe	236
हिंदी	252
عربي	270

MEET YOUR ZIIP HALO

Spark brilliant skin with ZIIP HALO - the only skincare device that works across an entire spectrum of electrical currents, allowing you to treat a whole range of skincare concerns, not just lifting and contouring.

Powered by Microcurrent and Nanocurrent™, each treatment has a unique combination of electrical waveforms to target a different area of the face or achieve a specific skin ambition.

WHAT IS IN THE BOX

ZIIP HALO
Silver Gel 30ml
Crystal Gel 30ml
Golden Gel 30ml
Charging USB-C cable
User Guide
Storage Bag

IMPORTANT SAFETY INFORMATION

- DO NOT operate the device where it could come into contact with water or steam.
- DO NOT attempt to take the device apart. The device is not serviceable by the user. This will void the warranty.
- STOP using the device immediately if you think it is defective.
- DO NOT use or charge the device with a power adapter that is damaged or faulty.
- DO NOT harshly bend or abuse the USB-C Charging Cord.

CUSTOMER CARE

Please contact ZIIP Beauty via the website if you have any questions or need any advice.

WARNINGS

ZIIP HALO is designed for cosmetic use only and for individuals in good health. If you have any medical concerns, such as a severe medical illness or have recently had surgery, please consult your doctor

before using ZIIP HALO.

Caution should be used for people with suspected or diagnosed heart problems - please consult your physician before use.

Caution should be used where there is a tendency to hemorrhage following acute trauma or on any areas of the skin which lack normal sensation.

Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium.

DO NOT hold stimulation over the temples.

ZIIP HALO should NOT be used during any activity in which involuntary muscle contractions may put the user at undue risk of injury.

ZIIP HALO should only be used on normal, intact skin. ZIIP HALO should NOT be used in the following areas:

- Breasts
- Eyelids
- Front of the throat (from under the chin to the collar bone)
- Genitals

If you experience irritation, extreme dryness, bumps or an allergic reaction, discontinue use immediately. Consult a physician if symptoms remain.

Electronic monitoring equipment (such as ECG monitors and ECG alarms) may not operate properly when a ZIIP HALO device is in use.

DO NOT submerge the device in water. If you store your ZIIP HALO in an area such as a bathroom, ensure the product does NOT come in contact with water.

DO NOT modify this device in any way. Doing so will invalidate the warranty and may cause harm. Do not use this device if it is damaged or has been submerged in water.

Keep out of reach of children.

CAUTIONS

Store the device at room temperature and out of direct sunlight.

Only use the charging cable supplied with the device.

Do not leave the device unattended while charging

Ensure the silicone port plug is inserted into the port while the device is being used

CONTRAINDICATIONS

ZIIP HALO should NOT be used by people who:

- are under the age of 18 years
- are pregnant
- are subject to seizures, including suspected or diagnosed
- epilepsy
- have cancer/tumors
- have a cardiac pacemaker with implanted defibrillators/ stimulators or have electronic implanted devices

The device should NOT be used on or near cancerous lesions, or applied over swollen, infected or inflamed skin e.g., phlebitis, thrombophlebitis, broken capillaries, varicose veins etc.

DO NOT apply stimulation across your chest as the introduction of electrical current into the chest may cause rhythm disturbances to your heart, which could be lethal.

DO NOT use this device on patients who have a cardiac pacemaker, implanted defibrillator, or other implanted metallic or electronic device, because this may cause electric shock.

SENSITIVITY TEST

1. Ensure your skin is clean and dry.
2. Apply a thumbnail amount of conductive gel to the treatment area and spread evenly.
3. Turn the device on and move in small circular motions for 3 minutes.
4. Stop the test immediately if you do not feel comfortable or if you see adverse effects.
5. Wait 48 hours and examine the tested areas. If the skin in the tested area appears normal, you can start a full treatment

DEVICE OVERVIEW

The ZIIP HALO device is an electrical beauty device that is intended for skin cell renewal and electrical muscle stimulation. It is indicated for at-home, over-the-counter cosmetic use.

HOW IT WORKS

The ZIIP HALO device delivers low-intensity electrical currents designed to mimic the body's natural bioelectrical signals. These currents are intended to support the skin's natural functions,

including the appearance of rejuvenated skin and improved facial tone. ZIIP HALO treatments use varying intensities of Microcurrent and Nanocurrent™ technologies. Clinically proven to show visible improvements in the appearance of fine lines and wrinkles, skin firmness as well as a reduction in visible blemishes, pigmentation and uneven skin tone, and increase skin radiance and glow.

- 1 Beep - On and ready to use
- 2 Beeps - Off
- 3 Beeps - Synced with app treatment
- Flashing Orange Light - Low Battery
- Orange Light (while plugged in) - Charging

HOW TO USE

CHARGING:

1. Charge the device before first use by connecting the charging cable to a USB power adapter and plugging it into your ZIIP HALO.
2. While charging, the band will glow orange. A full charge takes up to 4 hours and turns off when fully charged.
3. When the battery is low, the center band will flash orange.

CONNECT TO THE APP:

1. Download the ZIIP Beauty App from the Apple App Store or Google Play Store.
2. Open the app and select your preferred treatment.
3. Turn on your ZIIP HALO.
4. Tap “Start Treatment” in the app.
5. The device will beep and flash three times when synced successfully.
6. Follow the on-screen instructions.

INSTRUCTIONS FOR USE

PREPARE YOUR SKIN:

1. Wash and dry your face. Remove makeup, skincare products, and oil residue as this can prevent the device from working.
2. Take off jewelry from your face and hands.
3. Shave facial hair for better conductivity.

You must conduct a sensitivity test before first use.

1. Apply 5-6 pumps of any ZIIP Gel to clean, dry skin.
2. Single press the power button to turn on your device (it will beep once and glow).

3. The device will default to “The Lift” treatment. To access other ZIIP treatments, see app connectivity section above.
4. Place the globes onto your skin, holding the device so the charging port is facing down. The light band will brighten once the globes establish proper contact with your skin.
5. Starting at your jawline, glide the HALO outwards. The device will vibrate every 5 seconds to indicate when to move to the next section. When you hear a beep, switch to the other side of your face and repeat the process.
6. The treatment will run for 4 minutes and automatically turn off when complete.
7. Once finished, clean your device and wipe off any remaining gel from your face.
8. Place the device into the protective storage bag when not in use.

WHAT TO EXPECT DURING TREATMENT

During some treatments with higher intensities of Microcurrent, you will feel a tingling sensation. If this feels too intense for you, try speeding up your movement slightly or adding an extra pump of gel

to the treatment area.

It is normal to experience some mild redness on the skin during and just after each treatment.

FOR BEST RESULTS

For optimal results, use your ZIIP HALO 3–5 times per week. However, even once a week provides benefits.

Create a consistent treatment plan and follow it. If you miss a week, simply resume your routine to maintain progress

CLEANING, STORAGE & TRAVEL

Wipe the device with a soft lint-free cloth after use.

Avoid getting water, other liquids or excess conductive gel inside of any part of the device, including the USB charging port.

Do not use abrasive cleaner or chemicals on the device.

Keep the device stored at room temperature out of direct sunlight, and away from water, damp or debris that may damage the device.

Your ZIIP HALO is safe to travel with in the

United States, UK, EU and internationally. ZIIP accommodates 100-240VAC, 50/60Hz which is compatible with most countries in the world. The device comes with a standard USB cable; a country specific plug adapter is required for wall charging.

POSSIBLE ADVERSE EFFECTS

Mild adverse reactions include irritation and redness, but these should resolve quickly after a few hours.

Severe adverse effects are very rare, but if you experience any of the following during or after use, you should stop using the device immediately and consult your physician:

- severe skin irritation and burns
- persistent headache and/or other persistent painful sensations
- constant perception of flashing lights after you have finished treatment

DISPOSAL

The preferred method of disposal is through recycling. For environmental reasons, do not

dispose of the device in the household waste. Dispose of it at a suitable local collection or recycling point in accordance with EC Directive- WEEE (Waste Electrical and Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.

The box and user manual are card / paper and fully recyclable.

WARRANTY

Under normal use this product is guaranteed for 2 years after its original purchase date against manufacturing defects in materials or workmanship. This warranty extends only to the original purchaser. If the product should become faulty within the warranty period, contact The Beauty Tech Group.

- The Beauty Tech Group reserves the right to repair a defective product.
- The Beauty Tech Group reserves the right to replace a defective product with the most compatible product currently available.

This warranty does not cover damage caused by the following:

- Accident, misuse, abuse or alteration

- Servicing by unauthorised persons
- Connection to incorrect power sources
- Use with unauthorised accessories

The Beauty Tech Group shall not be responsible for any incidental, special or consequential damages resulting from the use of this product. All implied warranties, including but not limited to implied warranties of fitness and merchantability are limited in duration to two years from date of original purchase. This warranty gives you specific legal rights, and you may also have other rights that exclusion or limitation of incidental or consequential damages, or limitations on how long an implied warranty lasts, so some of the above limitations and exclusions may not apply to you. This warranty is given in lieu of all other warranties, written or oral, whether expressed by affirmation, promise, description, drawing, model or sample. Any and all warranties other than this one, whether expressed or implied, including implied warranties of merchantability and fitness for a particular purpose are hereby disclaimed. If the product is defective within the warranty period please contact The Beauty Tech Group.

DEVICE SPECIFICATIONS

Power levels:	900nA - 200mA
External power supply:	Input 100-240VAC, 50-60Hz, Output 9VDC, 1.5A
Electrical output:	3.7V
Device size and weight:	67.5g, 29mm x 50mm x 99mm 2.4oz, 3.91" x 1.98" x 1.15"

STORAGE INSTRUCTIONS

STORAGE AND TRANSPORTATION CONDITIONS

Temperature:	5°C to 45°C
Humidity:	0% and 95%
OPERATIONAL CONDITIONS	
Temperature:	5°C to 45°C
Humidity:	0 to 95%
Pressure:	3,000m.

SYMBOLS EXPLAINED

Model:	ZIIP HALO
Power:	9V, 1.5A, 100
Temperature:	5°C to 45°C

Use the device only with original power supply. For your safety, never use DC adapters which are not approved by the manufacturer.

Double insulation.

Consult instructions for use.

Waste electronic and electrical equipment.

Please read all warnings and safety precautions prior to use.

Protect from moisture.

This device must not be used near to water or in water itself (e.g., in a wash basin, in the shower or bathtub) - risk of electric shock!

Our products and packaging are recyclable, do not throw them away.

This symbol means: Conforms to certain EC Directives on health, safety and environmental protection standards. CE stands for 'Conformité Européenne'

This symbol implies conformity to UK Directives on health, safety, and environmental protection standards. UKCA stands for UK Conformity Assessed.

EMC INFORMATION

The following information is provided in accordance with EN

IEC 60601-1, IEC 60601-1-11, IEC 60601-2-57. ZIIP

HALO is compliant with IEC 60601-1-2

GUIDANCE AND MANUFACTURER'S DECLARATION ELECTROMAGNETIC EMISSIONS		
ZIIP HALO is intended for use in the electromagnetic environment specified below. The customer or user of ZIIP HALO should assure that it is used in such an environment		
Emissions	Compliance level	Electromagnetic Environment Guidance
RF Emissions CIS-PR 11	Group 1 Class B	ZIIP HALO uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment
Harmonic emissions IEC 61000-3-2		ZIIP HALO does not operate when connected to a power line
Voltage fluctuations / flicker emissions IEC 61000-3-3		ZIIP HALO does not operate when connected to a power line

GUIDANCE AND MANUFACTURER'S DECLARATION ELECTROMAGNETIC IMMUNITY			
ZIIP HALO is intended for use in the electromagnetic environment specified below. The customer or user of ZIIP HALO should assure that it is used in such an environment			
Immunity Test	IEC 60601 Test Level	Compliance level	Electromagnetic Environment Guidance
Electrostatic discharges (ESD) IEC 61000-4-2	±6 kV contact ±8kV air	±6 kV contact	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/burst IEC 61000-4-5	±2 kV for power supply lines ±1 kV for input/output lines	N/A	ZIIP HALO does not operate when connected to a power line
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	N/A	ZIIP HALO does not operate when connected to a power line
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 6100-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (95% dip in UT for 5 sec)	N/A	ZIIP HALO does not operate when connected to a power line
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment
NOTE: UT is the a.c. mains voltage prior to application of the test level			

GUIDANCE AND MANUFACTURER'S DECLARATION ELECTROMAGNETIC IMMUNITY			
ZIIP HALO is intended for use in the electromagnetic environment specified below. The customer or user of ZIIP HALO should assure that it is used in such an environment			
Immunity Test	IEC 60601 Test Level	Compliance level	Electromagnetic Environment Guidance
Conducted RF IEC 6100-4-6 Radiated RF IEC 6100-4-3	3 Vms 150 kHz to 80 MHz 3 Vms 80 MHz to 2.5 GHz	3 Vms 150 kHz to 80 MHz 3 Vms 80 MHz to 2.5 GHz	Portable and mobile RF communications equipment should be no closer to any part of ZIIP HALO, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter RECOMMENDED SEPARATION DISTANCE $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ (80 MHz to 800MHz) $d = 2.3 \sqrt{P}$ (800 MHz to 2.5 GHz) Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation in meters (m) Field strengths from fixed RF transmitters as determined by electromagnetic site surveys should be less than the compliance level in each frequency range. ¹⁾ Interference may occur in the vicinity of equipment marked with the following symbol 

NOTE 1: at 80 MHz and 800 MHz, the higher frequency range applies.
 NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/ cordless) telephones and land mobile radios, amateur radios, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy.
 To assess the electromagnetic environment due to fixed RF transmitters, and electromagnetic site survey should be considered. If the measured field strength in the location in which the ZIIP HALO is used exceeds the applicable RF compliance level above, ZIIP HALO should be observed to verify normal operations. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the ZIIP HALO.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.

RECOMMENDED SEPERATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE ZIIP HALO

ZIIP HALO is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of ZIIP HALO can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the ZIIP HALO as recommended below, according to the maximum output power of the communications equipment

Rated Maximum Output Power Of Transmitter W	Separation Distance According To Frequency Of Transmitter M		
	$d=1.2\sqrt{P}$	$d=1.2\sqrt{P}$	$d=1.2\sqrt{P}$
.01	.20	.20	.25
.1	.40	.40	.75
1	1.2	1.2	2.4
10	3.7	3.7	7.4
100	12	12	23.4

NOTE 1: at 80 MHz and 800 MHz, the higher frequency range applies.
 NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people

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ZIIP HALO is manufactured for **The Beauty Tech Group Ltd.**