# Electric Pad for Neck and Shoulders InnovaGoods -

## InnovaGoods®

### User's Manual

## ENGLISH

IMPORTANT INSTRUCTIONS. RETAIN FOR FUTURE USE.

### TECHNICAL SPECIFICATIONS

- DC In: 220-240 V / 50-60 Hz
- Max. power: 70 W
- 3 temperature levels
- Removable remote control
- Automatic switch off: 120 mins

### INSTRUCTIONS

- Plug the Electric Pad for Neck and Shoulders InnovaGoods® into a mains socket.
- Place it on your back and shoulders (see image). You can also use it on other areas of the body.
- Select the temperature by pressing the button on the remote control. The LED will indicate the temperature selected: «LO»: 40 °C, «MED»: 50 °C and «HI»: 60 °C (+/- 5 °C).
- If you need to wash the product, unplug the cable by pulling gently on the connector; open the zip, remove the cover, wash it and leave it to dry fully. The resistor inside MUST NOT be washed.

### WARNINGS

Do not use if wet. Do not use while sleeping. Do not leave the product unsupervised. Only use with the remote control included. Prolonged use could cause burns. Do not insert pins. Do not use folded or rucked. The appliance is an over-blanket. This appliance must not be used by persons insensitive to heat and other very vulnerable persons who are unable to react to overheating like children under the age of 3. The heating pad can be used by children older than 3 years and younger than 8 years under supervision and with the control always set to the minimum temperature value. For continuous use set the appliance to the minimum temperature value. When not in use, once at

room temperature, store the appliance in the packaging provided. Do not crease the pad by placing items on top of it during storage. This appliance is not intended for medical use in hospitals. Do not use if wet. Examine the appliance frequently for signs of wear or damage. Do not dry clean. Do not bleach. Do not iron. Do not tumble dry. The detachable cover can be washed by hand and shade dried. The cable and control unit of the appliance can give rise to risks of entanglement, strangulation, tripping or treading if not correctly arranged. Make sure that excess cable ties and electric cords are securely arranged. Make sure that the voltage of the power source corresponds to that of the appliance. Do not pull on the cable and ensure it does not tangle. Keep the cable away from hot and sharp surfaces. Place the appliance on a level, dry and heat-resistant surface. When you turn the device on for the first time, you may notice a strange smell. This is normal and temporary. Unplug the power cable when the appliance is not being used and let it cool completely before cleaning or handling. Clean with a damp cloth. Do not use abrasive products for cleaning. Do not immerse the appliance or cable in liquid. Do not handle it with wet hands Do not use near water or damp areas. Suitable for domestic use only. Use the device for its intended purpose only. This appliance is not a toy. Keep out of reach of children. Store the device in a dry, dust-free place. Do not leave the appliance unsupervised when it is turned on. Do not modify the appliance or any of its components. If the appliance, the cable or the pieces looked damaged or do not work correctly, do not use them or attempt to repair them; they should only be repaired by a qualified technician. This device is not intended for use by children or persons with reduced physical, sensory or mental capabilities, unless they are supervised by an adult responsible for their safety. This device is not household waste. Please follow the environmental regulations in place for its recycling Discard the product in accordance with the European Directive 2002/96/EC. Enquire about your local recycling legislation for electrical and electronic products and follow these steps. Do not dispose of this product with household waste. Following these instructions will help prevent harmful consequences to personal and environmental safety.

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