Arm and Leg Double Pedal Exerciser InnovaGoods -

ENGLISH

WARNINGS

IMPORTANT: KEEP THE INSTRUCTIONS FOR FUTURE REFERENCE. Use the appliance on a flat, stable surface. Check correct assembly before use. Base your training programme on your physical condition. Doing an inadequate warm-up or exercises which are too intense for your physical condition can cause damage or injuries. Do not mix new and used batteries in the compartment. If you are not going to use the unit for a long period of time, remove the battery. Do not expose batteries to high temperature or direct sunlight. Do not attempt to open them. do not dispose of them in fire, there is risk of explosion. Clean with a damp cloth. Do not use abrasive products for cleaning. Do not immerse the appliance in liquid. Do not handle it with wet hands Do not use near water or damp areas. Suitable for domestic use only. Use the device for its intended purpose only. Store the device in a dry, dust-free place. Do not leave the device unsupervised while in use. To avoid risks, do not use it if any parts are damaged and do not attempt to disassemble or repair it. This device is not intended for use by children or persons with reduced physical, sensory or mental capabilities, unless they are supervised by an adult responsible for their safety. This device is not household waste. Please follow the environmental regulations in place for its recycling Read the instructions carefully before first use. Keep the manual for further queries. Discard the product in accordance with the European Directive 2002/96/EC. Enquire about your local recycling legislation for electrical and electronic products and follow these steps. Do not dispose of this product with household waste. Following these instructions will help prevent harmful consequences to personal and environmental safety. This product contains a coin cell battery. Ingestion of a coin cell battery can cause internal chemical burns, soft tissue perforation, and death. Severe burns may occur within two hours of ingestion. Dispose of used batteries immediately. Keep new and used batteries out of reach of children. If you think the batteries have been swallowed or introduced to any part of the body, seek immediate medical attention. Do not allow children to change batteries without adult supervision.

CONTENTS

- 1 curved base (A)
- 2 leg supports (F,G)
- 2 poles (B,C)
- 2 pedal exercisers (D,E)
- 1 digital controller (LR44, 1.5 V battery)
- 10 nuts and screws (J)
- 2 assembly tools (H)

INSTRUCTIONS

- Follow the steps in the assembly guide to assemble the Arm and Leg Double Pedal Exerciser InnovaGoods®. Follow steps 1 to 8. Each image uses letters to show which pieces should be assembled using the nuts, screws (J) and the tools included (H).
- Once assembled, position the digital controller (9) on one of the pedal exercisers (D or E) to monitor the exercise of arms or legs.
- Adjust the height of the pole (B,C) and the closing pressure as shown in the image (5).
- Adjust the pedalling intensity of the

- hands or feet by adjusting the knobs (I).
- Digital controller functions: press the «MODE» button: 1-**CNT** (indicates turns of the handlebar; 2-**TIMER** (indicates the time of use); 3-**CAL** (indicates the calories consumed); 4-**DISTANCE** (indicates the distance in km: 0.00); 5-**TCNT** (indicates accumulated turns); 6-**SCAN** (shows each of the above values every 8 s).
- Press and hold the MODE button for 5 seconds to reset all of the values.
- If you need to change the battery of the controller, remove the front of it by pulling the lid.

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