ENGLISH

WARNINGS

IMPORTANT: KEEP THE INSTRUCTIONS FOR FUTURE REFERENCE. Do not wash this device in a washing machine or get it wet. Do not mix new and used batteries in the compartment. If you are not going to use the unit for a long period of time, remove the battery. Do not expose batteries to high temperature or direct sunlight. Do not attempt to open them. Do not dispose of them in fire, there is risk of explosion. Clean with a damp cloth. Do not use abrasive products for cleaning. Do not immerse the appliance in liquid. Do not handle it with wet hands. Do not use near water or damp areas. Suitable for domestic use only. Use the device for its intended purpose only. Store the device in a dry, dust-free place. Do not leave the device unsupervised while in use. To avoid risks, do not use it if any parts are damaged and do not attempt to disassemble or repair it. This device is not intended for use by children or persons with reduced physical, sensory or mental capabilities, unless they are supervised by an adult responsible for their safety. This device is not household waste. Please follow the environmental regulations in place for its recycling. Make sure that the voltage of the power source corresponds to that of the appliance. Unplug the appliance whenever you are not using it. Do not pull on the cable and ensure it does not tangle. Keep the cable away from hot and sharp surfaces. Before cleaning, make sure that appliance is disconnected from the power supply and has cooled completely. Clean with a damp cloth. Do not use abrasive products for cleaning. Do not immerse the appliance or cable in liquid. Do not handle it with wet hands. Do not use near water or damp areas. Suitable for domestic use only. Use the device for its intended purpose only. Store the device in a dry, dust-free place. Do not leave the appliance unsupervised when it is turned on. This appliance is not a toy. Keep out of reach of children. Do not modify the appliance or any of its components. If the appliance, the cable or the pieces looked damaged or do not work correctly, do not use them or attempt to repair them; they should only be repaired by a qualified technician. This device is not intended for use by children or persons with reduced physical, sensory or mental capabilities, unless they are supervised by an adult responsible for their safety. This device is not household waste. Please follow the environmental regulations in place for its recycling. Discard the product in accordance with the European Directive 2002/96/EC. Enquire about your local recycling legislation for electrical and electronic products and follow these steps. Do not dispose of this product with household waste. Following these instructions will help prevent harmful consequences to personal and environmental safety.

TECHNICAL SPECIFICATIONS

- Control module
- 4 inflatable leg wraps
- 1 set of connectable tubes
- AC converter: 110-240 V / 50-60 Hz
- DC In: 6 V / 800 mA
- Batteries: 4 x AA (1.5 V)
- Consumption: 6.5 W
- Pressure: 96 kPa
- 4 pressure intensities

INSTRUCTIONS

- Connect tube (1) to the control module. Connect tubes (2) to tube (1) and tubes (3) to tubes (2).
- Connect the ends of the tubes (2) to the upper leg wraps (2); connect the ends of the tubes (3) to the lower leg wraps (3) which each have 2 air inlets.
- Position the upper leg wraps (2) onto the thighs and secure them firmly with the velcro.
- Position the lower leg wraps (3) onto your legs and secure them firmly with the velcro.
- Connect the Air Compression Leg Massager InnovaGoods® to the mains electricity and the control module (4).
- Press the «ON/OFF» button to switch the device on.
- The Air Compression Leg Massager InnovaGoods® will start to inject air into the leg wraps; you should wait for a few minutes until they reach the selected pressure.
- It features 4 programmes from lower to higher intensity inflation/pressure. Press the «+» or «-» button to change the programme. Each programme has a different inflation/pressure time and relaxation time.
- Unplug the Air Compression Leg Massager InnovaGoods® from the mains when not in use.
- Using with batteries: remove the rear cover from the remote control and insert the batteries. Unplug the DC cable from the control module and switch the device on. NOTE: the Air Compression Leg Massager InnovaGoods® is recommended for use with the mains cable at all times, using it with batteries should only occur on occasions where it is not possible to plug the device into the mains electricity, as it will lose power if the batteries are not fully charged.

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