Eye Massager with Phototherapy InnovaGoods -

Clean the head with a towel or tissue paper after each use.

ENGLISH

WARNINGS

IMPORTANT: KEEP THE INSTRUCTIONS FOR

FUTURE REFERENCE. Fully recharge the battery the first time. Make sure that the voltage of the power supply is compatible with the device. Do not pull on the cord and make sure that it does not become tangled. Keep the cord away from hot or sharp surfaces. Clean the device with a slightly damp cloth. Do not use abrasive products or materials to clean it. Do not put the device or cord in liquid. Do not handle it with wet hands. Do not use the device close to a source of water or in damp places. Suitable for domestic use only. Use the device for its intended purpose only. Store the product in a dry place with no dust. If the device, the pieces or the cord show defects or they do not work properly, do not use them or attempt to repair them. Only a qualified technician should do it. This device is not designed for use by children or people with physical and mental disabilities, unless they are fully supervised by an adult who is responsible for their safety. This product is not household waste, so follow the environmental standards for recycling Discard the product in accordance with the European Directive 2002/96/EC. Enquire about your local recycling legislation for electrical and electronic products and follow these steps. Do not dispose of this product with household waste. Following these instructions will help prevent harmful consequences to personal and environmental safety.

TECHNICAL SPECIFICATIONS

- DC In: 5 V / 500 mA
- Battery: LiPo 3.7 V / 300 mAh
- Vibration: 1200 rpm
- Thermotherapy: heat (42 °C max.)
- Phototherapy: red and blue LED
- Touch vibration mode

INSTRUCTIONS

- Charge the Eye Massager with Phototherapy InnovaGoods® in a USB source. The red LED will flash while it is charging; once charged, the LED will turn blue and remain fixed.
- Press the middle button of the device for 2 seconds to switch on/off. Briefly press the

- button again to activate the modes: Press once: vibration; press twice: heat + red LED + vibration; press three times: blue LED + vibration.
- Make even movements from side to side around the area being treated.
- The vibration works automatically when the head of the massager and the base are in contact with the skin.
- Vibration mode: relaxes the eye muscles. It can be used along with cream/gel to look after the eyes, achieving a greater absorption of the product.
- Thermotherapy: the head will start to heat up until it reaches a maximum of 42 °C and will remain at a constant temperature.
- Phototherapy: Red LED for wrinkles and expression lines in the eye area. Blue LED for anti-fatigue and firming effect of the skin.
- The massager will automatically switch off after 2 minutes without use.

Clean the head with a towel or tissue paper after each use.

ENGLISH

WARNINGS

IMPORTANT: KEEP THE INSTRUCTIONS FOR

FUTURE REFERENCE. Fully recharge the battery the first time. Make sure that the voltage of the power supply is compatible with the device. Do not pull on the cord and make sure that it does not become tangled. Keep the cord away from hot or sharp surfaces. Clean the device with a slightly damp cloth. Do not use abrasive products or materials to clean it. Do not put the device or cord in liquid. Do not handle it with wet hands. Do not use the device close to a source of water or in damp places. Suitable for domestic use only. Use the device for its intended purpose only. Store the product in a dry place with no dust. If the device, the pieces or the cord show defects or they do not work properly, do not use them or attempt to repair them. Only a qualified technician should do it. This device is not designed for use by children or people with physical and mental disabilities, unless they are fully supervised by an adult who is responsible for their safety. This product is not household waste, so follow the environmental standards for recycling Discard the product in accordance with the European Directive 2002/96/EC. Enquire about your local recycling legislation for electrical and electronic products and follow these steps. Do not dispose of this product with household waste. Following these instructions will help prevent harmful consequences to personal and environmental safety.

TECHNICAL SPECIFICATIONS

- DC In: 5 V / 500 mA
- Battery: LiPo 3.7 V / 300 mAh
- Vibration: 1200 rpm
- Thermotherapy: heat (42 °C max.)
- Phototherapy: red and blue LED
- Touch vibration mode

INSTRUCTIONS

- Charge the Eye Massager with Phototherapy InnovaGoods® in a USB source. The red LED will flash while it is charging; once charged, the LED will turn blue and remain fixed.
- Press the middle button of the device for 2 seconds to switch on/off. Briefly press the

Clean the head with a towel or tissue paper after each use.

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