ENGLISH

TECHNICAL SPECIFICATIONS

- DC In: 5 V / 600 mA
- Battery: Li-ion 500 mAh
- Up to 6000 rpm
- 3 massage and heat modes
- Heat function
- USB cable

INSTRUCTIONS

- Charge the Double Chin Reducer with LED InnovaGoods® in a USB source; the LED on the head will flash while it is charging and will stop once charging is complete.
- Wash and moisturise the massage area with your usual cream. You can use it on the neck, face, arms, legs, etc.
- Press and hold the button (1) to switch the device on. Press the button (1) to change between the 3 massage modes.
- Hold and slide the metal sensor (2) across the treatment area. Make sure the hand holding the device is touching the metal sensor on the handle (3).
- Preventive massage (young skin): Blue LED + short intermittent vibrations + heat.
- Normal massage (middle-aged skin): Green LED + Long intermittent vibrations + heat.
- Repairing massage (mature skin or marked areas): Red LED + constant vibration + heat.
- Once finished, switch the device off (1) and clean it with a dry cloth.

WARNINGS

Do not use the product if you are suffering from any inflammation or lesions in the area. If the Double Chin Reducer with LED InnovaGoods® does not vibrate, make sure the hand holding the handle is touching the metal sensor (3).

Fully charge the battery before its first use. Make sure that the voltage of the power source corresponds to that of the appliance. Do not pull on the cable and ensure it does not tangle. Keep the cable away from hot and sharp surfaces. Clean with a damp cloth. Do not use abrasive products for cleaning. Do not immerse the appliance or cable in liquid. Do not handle it with wet hands. Do not use near water or damp areas. Suitable for domestic use only. Use the device for its intended purpose only. Store the device in a dry, dust-free place. Do not leave the appliance unsupervised when it is turned on. If the appliance, the cable or the pieces looked damaged or do not work correctly, do not use them or attempt to repair them; they should only be repaired by a qualified technician. This device is not intended for use by children or persons with reduced physical, sensory or mental capabilities, unless they are supervised by an adult responsible for their safety. This device is not household waste. Please follow the environmental regulations in place for its recycling. Discard the product in accordance with the European Directive 2002/96/EC. Enquire about your local recycling legislation for electrical and electronic products and follow these steps. Do not dispose of this product with household waste. Following these instructions will help prevent harmful consequences to personal and environmental safety.

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ALL RIGHTS RESERVED. MADE IN P.R.C.