MANUAL DE INSTRUCCIONES

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THERMOBI

Manual de instrucciones	ES
Instruction manual	EN
Manuel D'Utilisation	FR
Manuale D'Istruzioni	IT
Anweisungen	DE
Instruções	PT
Gebruiksaanwijzing	NL
Bruksanvisning	SW
Instrukcja obsługi	PL
Brugsanvisning	DA

ES - Lea las instrucciones antes de usar el producto.
EN - Please read instructions before use.
FR - Lisez les instructions avant d'utiliser le produit.
IT - Leggere le istruzioni prima di usare il prodotto.
DE - Lesen Sie die Gebrauchsanweisung, bevor Sie das Produkt verwenden.
PT - Leia as instruções antes de utilizar o produto.
NL - Lees de instructies voordat u het product gebruikt.
SW - Läs instruktionerna innan du använder produkten.
PL - Läs instruktionerna innan du använder produkten.
DA - Læs vejledningen, før du bruger produktet.

IMPORTANTE

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IMPORTANT

Keep the instructions for future reference. Read these instructions before assembling and using this product. Adult assembly is required.

Thank you for choosing the THERMOBI household ear thermometer! In this manual you will find all the information you need to use and maintain this product correctly and safely.

Before you start, we recommend that you read this manual carefully to ensure optimal and safe use of the product.

As a guarantee of the quality of our products, the manufacturer has a management system certified under ISO 13485. And this product complies with the requirements of Regulation (EU) 2017/745 on medical devices.

Notice to the user and/or patient: any serious incident that has occurred in relation to the device must be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is established.

1. WARNINGS AND PRECAUTIONS

1. WARNING! Do not disassemble, repair or modify this device. It may cause errors in the measurement results or failure of the device.

2. WARNING! The infrared ear thermometer is an accurate device; do not allow children to play with it. Be careful not to drop it from a height, as a fall may loosen the sensor and affect the accuracy of the temperature measurement.

3. WARNING! The infrared ear thermometer must not be transported without its packaging to prevent mechanical vibrations from affecting the accuracy of the temperature measurement.

4. WARNING! If it is not used for a long time, remove the batteries and keep it out of the reach of children. Also, if the battery is not used for a long time, it may be damaged due to leakage.

5. WARNING! This product contains sensitive electronic components and should be avoided in environments with strong electromagnetic interference (such as mobile phones, microwaves, etc.) to prevent temporary impacts on its accuracy.

6. WARNING! Read this instruction manual carefully before using the device and confirm that the batteries are installed.

7. WARNING! Keep the probe clean before and after use. When the sensor lens becomes dirty, gently wipe it with a dry cloth or cotton swab. Do not wipe it with other objects or blow on the infrared sensor with your mouth as this may scratch the sensor mirror or cause the device to malfunction.

8. WARNING! The materials of the application components of the device are certified to be biocompatible.

9. WARNING! The patient may use the thermometer to measure themselves, or someone else may use it to measure the patient. In either case, the methods of use and maintenance are the same.

10. WARNING! All service and maintenance must be carried out before or after use.

11. WARNING! The thermometer must not be serviced or maintained while in use.

PRECAUTIONS

The temperature may vary slightly from the left ear to the right ear due to various factors. It is therefore best to always measure the same ear for comparison.

- Measurement results are for reference only and cannot replace a doctor's diagnosis. Some people do not have a fever when they are ill. If you feel unwell, regardless of the measurement results, it is recommended that you consult a doctor immediately.
- Sleeping on the side can cause the temperature of the pressed ear to rise. It is best to wait a few minutes before measuring.
- The ear canal must be kept clean, as excess ear wax will affect accuracy.
- Do not attempt to measure when the device is wet, as this may give inaccurate results.
- Do not measure body temperature within 30 minutes after exercise, bathing or eating. At these times the body temperature deviates from normal.
- If there is a temperature difference between the storage location of the device and the measurement environment, the device must be placed in the environment of use for more than 30 minutes, otherwise the measurement results may be in error.
- Do not use the thermometer if there is blood or drainage in the external ear canal.
- If ear drops or other medications have been placed in the ear canal, use the ear without the medication.
- Complete occlusion of the ear canal due to cerumen (wax) may result in lower temperature readings.
- Do not leave the thermometer unattended near children. Small objects may present a choking hazard.
- If multiple measurements are required for the same person, remove the thermometer from the ear canal after each measurement and follow the measurement steps for retesting.
- Limitation: The device cannot replace a rectal thermometer for measuring rectal temperature.
- When measuring, make sure that the probe is aligned with the eardrum and in close contact with the ear canal.

2. PURPOSE OF THE PRODUCT

The intended purpose of the device is to display body temperature measured in the ear cavity by thermal radiation.

3. PRODUCT FEATURES


The ear thermometer is equipped with an infrared sensor, which can transform the infrared light emitted by the human ear membrane into a corresponding electrical signal. The signal is corrected by the amplifier and signal processing circuit according to the internal algorithm of the instrument and the target emissivity, and then transformed into the measured human body temperature value.

PRODUCT NAME	Household ear thermometer
MODEL	THERMOBI
MODEL REFERENCE	X-00001/65
PRODUCT WEIGHT	0,13 lb
PRODUCT DIMENSIONS	5,70"*1,45"*2,24"
MATERIALS	Polypropylene (PP)


OTHER CHARACTERISTICS

- Power supply: DC 3 V (2 x AAA 1.5 V alkaline batteries). Not included.
- Measuring range: 34.0°C~42.2°C (93.2°F~108.0°F).
- Temperature units: °C/°F.
- Accuracy: Within the temperature display range of 35.0°C - 42.0°C, (95.0°F - 107.6°F), the accuracy is $\pm 0.2^{\circ}\text{C}$ ($\pm 0.4^{\circ}\text{F}$). Outside the temperature display range of 35.0°C - 42.0°C (95.0°F - 107.6°F) the thermometer will display Lo or Hi, when lower or higher than the visible measurement range.
- Time interval between measurements ≤ 12 seconds.
- Automatic switch-off time: 60s \pm 10s.
- Measurement duration ≤ 2 seconds.
- Applied part type BF (corresponding to applied parts that are in contact with the patient for short or long periods of time): The probe.
- Degree of protection provided by the enclosure (IP code): IP22.
- Mode of operation: continuous operation.
- Device type: Non AP/APG (cannot be used in the presence of flammable anaesthetic gas mixed with air, oxygen or nitrous oxide).
- Lifetime: 5 years (excluding vulnerable parts and consumables).

4. PARTS AND COMPONENTS

DISPLAY

5. INDICATIONS AND CONTRAINDICATIONS

INDICATIONS

It is intended for use in infants (except premature or small for gestational age infants), children (under adult observation) and adults.


CONTRAINDICATIONS

Do not use the ear thermometer if there is inflammation in the ear canal. It is suggested to use the ear thermometer after recovery from surgery or trauma to the ear canal.

6. INSTRUCTIONS FOR USE

Battery installation:

The device is powered by 2 AAA alkaline batteries (not included). Lift the tab down and remove the battery cover, place the batteries in the compartment following the polarity indications (negative and positive) inside. Close the battery cover (the product will perform a functional self-test of the system, such as full display).

Product use:

1. Remove the cover and press the measuring button "  ", the backlight of the product will be fully displayed, as in figure (1), and it will enter the standby state for measurement after 1 second, as in figure (2).


Figure 1

Figure 2



2. During the measurement, turn the head to the side, leaving the ear upwards, and gently pull the ear backwards to straighten the ear canal. In children older than 1 year and adults, pull the ear backwards and upwards. Gently insert the probe of the infrared ear thermometer into the ear canal, ensuring a tight seal to measure the correct temperature of the eardrum.

3. Press and release the measurement button "  ". After approximately 1 second, the infrared ear thermometer will beep indicating the end of the measurement and display the results.

Note: If the system is busy processing data, it will display "  ". When "  " is displayed, you can take the next measurement.  

FUNCTIONS

1. High temperature alert: If the measured temperature exceeds 37.5°C, a long beep followed by three short beeps will sound to alert the user.
2. Memory function: Press the memory button "  ", the memory value can be queried. The infrared ear thermometer can store 7 memory values. If this specified number is exceeded, the most recent value will overwrite the oldest value.
3. In memory mode, press and hold the memory button "  " for 5 seconds, when CLr appears on the display, the memory will be cleared.
4. Auto power off function: The device will automatically turn off after 60 ±10 seconds of non-operation.
5. Manual power off function: Press and hold the measurement button "  " to turn it off.
6. Battery replacement: if only the "  " symbol is shown on the display, replace the batteries as soon as possible, if the "  " symbol is shown, replace the batteries immediately.
7. Temperature unit change: In off state, press and hold the measuring button "  " for 8-12 seconds to enter the unit change mode. Briefly press the memory button "  " to select °C or °F, then press the measurement button "  " to confirm.

7. MAINTENANCE AND CLEANING

MEASUREMENT, TRANSPORT AND STORAGE ENVIRONMENT

1. Measurement environment:

- Ambient temperature: +10°C ~ +40°C (50°F ~ 104°F).
- Relative humidity: 15 % to 90 %, non-condensing.
- Atmospheric pressure: 70 kPa ~ 106 kPa.

2. Transport and storage environment:

- Ambient temperature: -20°C ~ +55°C (-4°F ~ 131°F).
- Relative humidity: 15 % to 90 %, non-condensing.
- Atmospheric pressure: 70 kPa ~ 106 kPa.

CARE AND CLEANING

- Before each use, check the sensor lens for dirt. If it is dirty, gently wipe it with a cotton swab or soft cloth dampened with clean water. Make sure the probe is completely dry before the next measurement.
- Use a clean, soft cloth dampened with clean water to gently wipe the LCD screen and casing. Warning: Do not immerse this product in medical alcohol or other liquids.
- Clean devices should be stored in a dust-free and dry place.
- Avoid direct sunlight. Do not store the device in a place with high temperature, humidity, dust or corrosive gases. Aging of the sensor may affect the measurement accuracy.

BATTERY REPLACEMENT

Battery replacement: When the " " symbol is displayed, replace with 2 new AAA alkaline batteries, slide open the battery cover and remove the old batteries, insert the new batteries, making sure to align them correctly according to the indications inside the battery compartment.


If the product will not be used for an extended period of time, remove the batteries to prevent damage from battery leakage.



To protect the environment, please dispose of used batteries at appropriate collection points according to national or local regulations.

This product contains recyclable batteries and electronic waste. Do not dispose of in household waste, but take it to appropriate local collection points.

8. SYMBOLS RELATING TO PRODUCT SAFETY REQUIREMENTS AND THEIR MEANINGS

SYMBOLS	IMPLICATION
	Type BF applied part
	Caution
	Recyclable
	Manufacturer
	Authorized representative in the European Community / European Union
	Refer to instructions manual
	Symbol for the marking of electrical and electronic devices according to directive 2012/19/EU
	Safety and environmental protection use period for 10 years
IP22	Protection from ingress of particulates than $\geq 0,49''$. Dripping water falliing within 15° of vertical will not have a harmful effect on the infrared ear thermometer per IEC 60529
	Temperature limit

	Humidity limitation
	Atmospheric pressure limitation
	This way up
	Fragile, handle with care
	Keep dry
	This device fulfils the provisions of regulation (EU) 2017/745 (Medical Device Regulation)
	Serial number
	Date of manufacture
	Medical device
	Importer
	Batch code
	Unique device identifier
	MR unsafe

9. CONDITIONS OF GUARANTEE

This product is covered by a 3-year warranty against manufacturing defects.

Thank you for trusting Mobiclinic with our THERMOBI ear thermometer for home use! We value your trust in our products and are committed to always offering you the best quality.

If you have any questions or need further assistance, please do not hesitate to contact us through our customer service at clientes@mobiclinic.com. We appreciate your comments and suggestions for further improvement.

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CO., LTD.
GRUPO R. QUERALTO S.A.