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EN- Please read these safety instructions carefully before use, and keep them for future question.

1. Safety indications

- ♥ ALWAYS consult with your physician or supplier to determine proper adjustment and use of the device.
- ♥ Don't store outside and don't let under the sun (Beware : don't let the product in a car in case of high outside temperature !!)
- ♥ Don't use outside but only indoor.
- ♥ Check carefully all parts of this item (screws, tubes, clips..) and if all screws are tightened.
- ♥ To ensure that the safety level is kept to the highest possible standard, determined by its construction, this product should be serviced regularly (once a year) by specialist retailers.
- ♥ Before beginning your program of exercise, consult your doctor to ensure that you are fit enough to use the equipment. Base your program of exercise on the advice given by your doctor. Incorrect or excessive exercise may damage your health!
- ♥ Any interference with parts of the product that are not described within the manual may cause damage, or endanger the person using this machine.
- ♥ In choosing the location of the apparatus, ensure a sufficient safety distance from any obstacles. The apparatus must not be mounted in the immediate vicinity of main passageways (paths, doorways, corridors).
- ♥ The machine is designed for use by adults and children should not be allowed to play with it.
- ♥ Use only as a pedal exerciser and in proper way.
- ♥ No modification must be made on this device. Use only original parts to replace the broken or defect one and let this be done by the technical staff.
- ♥ After assembling the product, make sure that it is securely locked and fixed before using.
- ♥ In case you didn't use your product for a long time, let it checked by the technical staff, before using.
- ♥ Always check stability before using.
- ♥ Follow the advice of the professional who prescribed this medical device, respect duration and frequency of use.
- ♥ In case of discomfort consult your professional.
- ♥ Warm up 5 minutes before using the training device.
- ♥ Keep fingers, clothes, jewelry or items away from moving parts of the pedal exerciser (rotating pedals).

2. Intended use, indications and contraindications, installation

Intended use:

The pedal exerciser allows the maintenance and increase of muscle strength and mobility. It also supports the physiotherapeutic treatment and active care of bedridden patients. It can be used with the arms or legs, at home, in the office or in nursing homes and hospitals.

Indication:

A Pedal exerciser can be used for upper and lower body exercises to improve circulation, muscle strength, joint range of motion, and coordination in the home, in the office, or in the clinic.

Contraindications:

Heart problems - knee joint pathologies

Only specialists can adjust the product and teach you how to use it. A wrong use or adjustment could lead to overstress of the device or of the patient. The device must be adjusted individually by the technical staff to suit your needs. Ask the technical staff how to use this product and make sure that you understood clearly their explanation. In case of doubt, don't hesitate to ask again.

1-Installation & Description:

We propose two versions with and without display. Both version fold easily for ease of storage and transportation, feature four anti-slip rubber pads **B** on the feet to prevent sliding and to protect surfaces and the pedals come complete with adjustable toe straps **A**.

The version with display shows speed, distance, total time, and calories burned.

1. Unfold the 4 legs outwards until the locking pin clicks automatically in the hole, follow photo [A].
2. Turn the tension knob to adjust the paddle resistance. If you turn in the direction of clockwise, the resistance will be higher, follow photo [B].
3. Place the product on a table or a bed when training the upper body and on the floor when training the lower body.

Use:

Warning: the user and the device must be in a secured, stable position

when using the device. Slipping risks!

1-Adjust the toe straps **A** and start pedaling.

Digital Display (depends of the chosen model)

1 RPM	2 SCAN	3 CNT	4 TIME	5 CAL	6 MODE /RESET
revolutions per minute	scan the various functions	revolutions count	time in use	calories used	start-Operating modes

SCAN : Rotates through the indicated displays pausing for 5s on each. To reset the display press and hold the red button "MODE/RESET" **6** for 4s till display shows all zero's.

After a inactivity of 4s the STOP indicator will appear and the timer will stop. The "Automatic Power-Off" feature will turn the product completely off after a 2 minutes inactivity.

Use with the display:

1. To start using the product, first remove the element **7** and then press and hold the red button "MODE/RESET" **6** for 4s till display shows all zero's.
2. Press and hold the red button "MODE/RESET" **6** to select the desired function. Start pedaling.

Battery replacement

If the digital display begins to fade you will need to replace the battery. To do so follow the instructions below.

1. Carefully remove the display from the body of the exerciser.
2. The battery holder is located at the back of the display. Slide the display panel **8** downwards and remove the battery from its holder.
3. Replace with a new 1.5V AG13 battery and slide back the display panel.
4. Mount back the display on the pedal exerciser.

3. Cleaning

This product has to be cleaned regularly with a mild cleaning agent and a soft cloth. It may be washed with soap and water or disinfected with non-chlorinated products. The product must be dried after each cleaning.

4. Materials

The product is made out of chromed steel and various plastic parts.

5. Disposal

Dispose of this product in accordance with your regional waste disposal regulations.

6. Technical specifications

This product has a 2 years warranty, from the date of purchase, against any manufacturing defect. Duration of use : 2 years.

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